QBC[®] Malaria Test

4277-112-350

Rev C

English:	pages	1 – 6	Deutsch:	Seiten	18 – 24
Français:	pages	6 – 12	Italiano:	pagine	24 – 29
Español:	páginas	12 – 18	Portuguese:	pages	30 – 35

See symbol glossary at end of insert. / Voir le glossaire des symboles à la fin de la notice. / Siehe Symbol-Erklärungen am Ende der Packungsbeilage. / Vedere il glossario dei simboli alla fine del foglio illustrativo. / Consulte o glossário de símbolos no fim do folheto informativo. / Consulte el glosario de símbolos al final del prospecto.

INTENDED USE

The **QBC[®]** Malaria Test is a qualitative screening method for rapidly detecting the presence of malaria in centrifuged capillary or venous blood.

SUMMARY AND EXPLANATION


Throughout much of the world, the diagnosis of many blood parasites is generally established by microscopic examination of stained blood films. Direct observation of centrifuged blood in specially coated **QBC** tubes now permits field detection of malaria parasites in a 10 sec to 3 min examination of parasite-infected blood. Parasite diagnosis by the **QBC** method is more sensitive than conventional blood-film microscopy due, in part, to the comparatively larger volume of blood screened for parasitic infection.¹ The **QBC** method additionally concentrates most of the parasites into the narrow zone of the blood tubes that is optically examined. Components of the **QBC** Malaria Test require no controlled test environment or refrigeration facilities. Centrifuged **QBC** Malaria Test tubes are examined under a fluorescence microscope, or a light microscope using a **QBC ParaLens Advance[™]** LED fluorescence microscope attachment.

PRINCIPLES OF THE PROCEDURE

QBC centrifugal hematology analysis utilizes density gradient layering of stained blood cells, together with mechanical expansion of the microhematocrit buffy coat. A natural extension of this **QBC** technology is provided by the less dense, parasite-infected cells concentrated within the separated red cell and white cell layers of the spun **QBC** tube. **QBC** Malaria Test tubes are internally coated with acridine orange stain and anticoagulants. The use of acridine orange dye to visualize malaria parasites has been shown to have advantages.¹⁰ Expansion of the centrifugally separated cell layers is achieved with a plastic float, inserted into the tube after filling with capillary or venous blood and sealing with a plastic closure. Separated, expanded layers, clearly identifiable in the spun tube, are dark and light red cells, granulocytes, lymphocytes-monocytes, and platelets, as illustrated below. Examination under a fluorescence microscope of the interface region between the light red cells and granulocytes (where the parasites are characteristically most abundant^{2,3}) permits rapid and unambiguous detection of the distinctively stained malaria organisms.

REAGENTS

Each **QBC** Malaria Test tube is internally coated with Acridine Orange, Potassium Oxalate, Sodium Heparin and K₂EDTA.

WARNINGS AND PRECAUTIONS

For *in vitro* Diagnostic Use.

- Malaria Tubes can be used with capillary blood collected from the finger or heel, or with venous blood.
- Always wear protective laboratory gloves when collecting and handling blood.

CAUTION

BLOOD SPECIMENS MAY CONTAIN THE HUMAN IMMUNODEFICIENCY VIRUS (HIV), HEPATITIS B VIRUS OR OTHER INFECTIOUS AGENTS. HANDLE BLOOD AS A POTENTIAL BIOHAZARD, CAPABLE OF TRANSMITTING INFECTION.

STABILITY AND STORAGE OF QBC TUBES

QBC Tubes may be used for malaria testing up to the expiration date stamped on the package label. Store the Tubes at 16°C to 37°C (60° to 99°F), away from heat and direct sunlight. Do NOT expose tubes to sunlight or overhead light for extended periods since the coating of acridine orange stains may deteriorate.

BLOOD COLLECTION

For malaria diagnosis, use blood from a finger or heel puncture, or a freshly-collected venous specimen. Observe the following general precautions:

- Be sure to clean the skin area thoroughly with antiseptic agent and dry with a sterile wiper.
- For capillary blood, puncture the skin with a sterile lancet. Collect 55 to 65 μL of blood directly into the **QBC** Tube as described under TEST PROCEDURES.
- For venous blood, draw the specimen with a sterile syringe or evacuated collection device. See TEST PROCEDURES for filling the **QBC** Tube.

CAUTION

TO PREVENT TRANSMITTING INFECTION, BE SURE THAT ALL BLOOD COLLECTION MATERIALS ARE THOROUGHLY STERILIZED BEFORE USE.

Anticoagulants: **QBC** Malaria Tubes are internally pre-coated with sodium heparin and dipotassium EDTA.

Interfering Substances: Leukocyte nuclei or the presence of multiple artifacts in blood do not appear to confuse or mask the detection of parasites by the **QBC** method.

Stability of Blood-Filled Tubes: Blood-filled **QBC** Malaria Tubes should be centrifuged promptly after preparation (see TEST PROCEDURES). After centrifugation, tubes can be held for as long as 3 days unrefrigerated, at temperatures up to 37°C (99°F), before microscopic examination. During storage, tubes must be kept upright, with the

closure end down. Tubes refrigerated at 4°C (39°F) can be held for up to two weeks prior to examination.

TEST PROCEDURES

Materials Provided: Malaria Tubes, Floats, Closures, Package Insert, Pressure-Sensitive Tube Labels for patient identification, QBC holographic stickers.

Materials Required But Not Provided:

- Sterile Lancets
- Forceps
- Sterile Wipes
- Fluorescence Optical Immersion Oil
- **QBC Capillary Centrifuge**
- **QBC Work Station**
- **ParaViewer** tube holder, and
- **ParaLens Advance™** LED fluorescence microscope attachment, with 60x objective, and compound light microscope with 10x High Point, Wide-Field Eye Piece(s), or
- Standard Fluorescence Microscope, equipped with FITC filter cluster and 50x or 60x oil objective lens with working distance greater than 0.32 mm.

BLOOD TUBE PREPARATION AND TEST PROCEDURES


Each **QBC** test requires 55 µL to 65 µL of blood from a finger or heel puncture or from a freshly collected venous sample. Tube preparation and test procedures are described below.

CAUTION: Wear protective gloves when collecting and handling blood to avoid the transmission of disease.


PREPARE AND CENTRIFUGE BLOOD TUBE

1. Fill the **QBC** Malaria Tube, from end nearest the two blue lines, directly from a finger (or heel) puncture or a collection tube of well-mixed venous blood. Fill the tube by capillary action to a level between the two blue lines.


With lint-free tissue, wipe off any blood on the outside of the tube.


2. KEEP THE TUBE NEARLY HORIZONTAL and roll between the fingers several times to mix the blood with the anticoagulant coating.


3. TURN THE TUBE AROUND AND TILT, allowing the blood to flow to the end with the orange-coated stain. Roll the tube between the fingers 5 times to mix the blood with the staining agent.

4. TILT THE TUBE SLIGHTLY so that blood flows away from the orange-coated end by at least 1/4" (6 mm) to allow space for installing the closure; then place the index finger over the end of the tube nearest the blue fill lines.


5. PRESS A PLASTIC CLOSURE onto the unsealed end of the tube. Then manually twist and press the closure to form a tight seal.

CAUTION

DO NOT FORCE THE CLOSURE ONTO THE TUBE. TUBES ARE MADE OF GLASS AND MAY BREAK.

6. WITH A CLEAN FORCEPS, PICK UP A FLOAT and insert it into the unsealed end of the tube. Tap with the forceps until the float is inside the tube.

7. For identification purposes, record the patient number on a pressure-sensitive label supplied in the Kit. Then wrap the short end of the label tightly around the tube, being careful to position the label **between the two white lines on the tube**.
- Note:** Prepare other tubes for centrifugation, as required. Up to **20 tubes** can be centrifuged at one time.
8. PLACE THE TUBES into slots of the centrifuge rotor. TUBES MUST BE PLACED SO THAT THE ROTOR IS BALANCED. Install the rotor cover: then close and latch the lid. Start the centrifuge.
9. AFTER 5 MIN CENTRIFUGATION and when the rotor has stopped, remove the tubes from the centrifuge. Store the tubes vertically, closure end down, in the rack.

IMPORTANT: After centrifugation, QBC Malaria Tubes can be stored up to 3 days at 16°C to 37°C (60°F to 99°F) before microscopic examination for parasites; or tubes can be refrigerated at 4°C (39°F) for up to two weeks before examination.

MICROSCOPIC EXAMINATION FOR MALARIA

Using a fluorescence microscope or a light microscope fitted with a **ParaLens Advance** and a **ParaViewer** tube holder:

1. Insert the centrifuged **QBC** Tube into the groove of the **ParaViewer**. Position the Tube so that the sealed end extends over the depressed area of the holder.
2. Place the **ParaViewer** with **QBC** Tube on the stage of the microscope. Use a 10x, wide-field, high-point eyepiece(s) and a 50x to 60x objective lens with a minimum working distance of 0.34 mm. Add 2 – 3 drops of fluorescence optical immersion oil. Bring the buffy coat area of the Tube into focus.
3. Move the microscope stage until the top of the expanded RBC layer (near the RBC-Granulocyte interface) is in the field of view.
4. Examine the entire circumference of the tube by rotating the tube as it rests in the **ParaViewer**.
5. Examine the area from the RBC-Granulocyte interface towards the closure to a distance of 10 to 15 fields.
6. Look for schizonts and gametocytes of *Plasmodia* at the Granulocyte- Lymphocyte/ Monocyte interface for a maximum of 10 fields.
7. Total examination time to exclude a negative is approximately 2 min.

The presence of malaria parasites is indicated by the distinct bi-colored signet ring forms of trophozoites strikingly apparent in cells near the granulocyte layer. A central area, where the “stone” of the signet ring would be located will be **green**; the ring portion will be **green to red**.

Gametocytes of *P. falciparum* will appear as yellow sickle-shaped bodies. Schizonts of *P. vivax* can be recognized by the presence of malaria pigment which appears as a dark brown color.

Developing merozoites in mature schizonts will often be seen very clearly surrounding the dark brown pigment.

QBC HOLOGRAPHIC STICKER

Attach **QBC** holographic sticker to all reports of the malaria tests performed using the **QBC** Malaria Test. Note: The use of the **QBC** holographic sticker ensures the test reviewer that the test was performed using authentic **QBC** technology.

TEST LIMITATIONS

The **QBC** Malaria Test is a qualitative method for diagnosing “positive” or “negative” patients. While numerous reports have demonstrated that the **QBC** test is significantly more sensitive for detecting malaria than stained blood films,^{1,4-8} **QBC** may not detect very low levels of parasitemia. Experienced microscopists can discriminate between various species of malaria parasites using the **QBC** test.^{1,4,9} Examination of a thin blood smear may be required to differentiate species in some cases.

PERFORMANCE CHARACTERISTICS

The **QBC** (quantitative buffy coat) technique was compared to the Giemsa-stained thick blood film (GTF) under field conditions. Three hundred sixty-four (364) blood samples were collected and each was tested with both **QBC** and GTF techniques. The GTF technique indicated 86 positive blood samples and 278 negative samples. The **QBC** technique indicated 89 positive samples and 275 negative samples. Relative to the results obtained with GTF, the **QBC** technique had a sensitivity of 87.2% and a specificity of 95.0%; concordance between the tests was 93.1%.

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REFERENCES

1. Spielman, A., *et al.*, 1988: Malaria Diagnosis by Direct Observations of Centrifuged Samples of Blood. *Am. J. of Trop. Med. And Hyg.*, 39 (4), 337-442.
2. Worth, R.M., 1964: The Heparinized Capillary Tube as an Epidemiological Tool. LI. Concentration of Blood Parasites by Centrifugation. *Am. J. Hyg.*, 80,70-74.
3. Woo, P.T. and Hauck, L. 1987: The Hematocrit Centrifuge Smear Technique for the detection of mammalian Plasmodium. *Trans. Roy. Soc. Trop. Med. Hyg.*, 81, 727-728.
4. Rickman, L., *et al.*, 1989: Rapid Diagnosis of Malaria by Acridine Orange Staining of Centrifuged Parasites. *Lancet* (1989) , 68-71.
5. Levine, R.A., *et al.*: Detection of haemoparasites using quantitative buffy coat analysis tubes. *Parasit Today*: 1989; 5 (4), 132-134.
6. Parzy, D., *et al.*: Quantitative buffy coat test (**QBC** test), Monofluo Kit Falciparum, comparative value in the rapid diagnosis of malaria. *Med. Trop.*: 1990, 50:97-101.
7. Namisiripongpun, V., *et al.*: The acridine orange stained capillary tube (the **QBC** system) in diagnosis of malaria: a field trial. Report presented at the 2nd Western Pacific Congress on Infectious Diseases and Chemotherapy. December 11-14, 1990. Jomtien Pattaya, Thailand.
8. Pornsilaptip, J., *et al.*: Detection of Plasmodia in acridine orange stained capillary tubes (the **QBC** system). *Asian J. Trop. Med. Pub. Hlth.*: 1990, 21 (4), 534540.
9. Atlas of Blood Parasite Morphologies: Becton Dickinson and Company: 1990.
10. Shute, G., Sodeman, T.: Identification of Malaria Parasites by Fluorescence Microscopy and Acridine Orange Staining. *Bull. Wld. Hlth. Org.*: 1973, 48, 591-596.
11. Wang, X., S. Zhu, Q. Liu, A. Hu, Z. Zan, Q. Yu, and Q. Yin. 1996. Field evaluation of the **QBC** technique for rapid diagnosis of vivax malaria. *Bull. World Health Organ.* 74 (6):599-603.

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Prinsessegracht 20
2514 AP The Hague
The Netherlands

QBC Europe

A division of Woodley Equipment Company Ltd
Old Station Park Buildings, St Johns Street, Horwich, Bolton,
Lancashire, BL6 7NY, United Kingdom
01204 669033 and 01204 460446,

QBC Logo, QBC, ParaLens, ParaLens
Advance and ParaViewer are
trademarks of Drucker Diagnostics ©
Drucker Diagnostics

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Philipsburg, PA 16866 U.S.A.