QBC™ AccuTube

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Italiano: pagine 6 – 8
Español: páginas 8 – 10
Portuguese: pages 10 – 11

INTENDED USE

QBC AccuTubes are used with the **QBC Autoread Plus** System to provide a diagnostic hematology profile. Venous or capillary blood can be tested. The following quantitative values are obtained from a single AccuTube: hematocrit, hemoglobin, mean corpuscular hemoglobin concentration (MCHC), platelet count, white blood cell count, granulocyte count (% and number), and lymphocyte/monocyte count (% and number).

PRINCIPLES OF THE TEST

When a tube of blood is placed in a centrifuge and spun, the different types of cells in the blood separate into layers from the heaviest to the lightest. A special coating in the AccuTube stains the different cell layers. These layers are made larger by a plastic float inside the AccuTube. The **QBC** analyzer measures the cell layers and converts the measurements into reported values.^{1,2,3}

SUMMARY AND EXPLANATION

The **QBC Autoread Plus** System includes an analyzer, power pack, printer, centrifuge and test accessories. The system uses specially designed tubes known as the **QBC AccuTube**.

Each AccuTube has two black fill lines, a white identification line, reagents, and a stopper. The tube is filled by capillary action or with a **QBC** Pipetter, modified with a spacer, to aspirate 70 µL of venous blood. After mixing the tube and seating the stopper, a plastic float is inserted in the AccuTube. The tube is spun in a centrifuge for 5 minutes to separate the blood components. The tube is then placed in the **QBC Plus** Analyzer which displays and prints results.

REAGENTS

Each AccuTube contains a maximum of 6.5µg Acridine Orange, 0.65 mg Potassium Oxalate, 1.5 USP units Sodium Heparin (anticoagulant), 0.33 mg Potassium EDTA (anticoagulant), agglutinating agents and stabilizers.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

WARNING: ACRIDINE ORANGE REAGENT MAY BE TOXIC; DO NOT INGEST. AVOID CONTACT WITH SKIN, EYES, AND CLOTHING.

WARNING: BLOOD SPECIMENS MAY CONTAIN THE HEPATITIS B VIRUS (HBV), HEPATITIS C VIRUS (HCV), HUMAN IMMUNODEFICIENCY VIRUS (HIV), OR OTHER DISEASE CAUSING AGENTS. HANDLE ALL PATIENT SPECIMENS AS POTENTIAL BIOHAZARDS CAPABLE OF TRANSMITTING INFECTION. WEAR APPROPRIATE PERSONAL PROTECTIVE EQUIPMENT, INCLUDING GLOVES, WHEN COLLECTING AND PROCESSING BLOOD.

WARNING: **QBC** ACCUTUBES ARE MADE OF GLASS. BE CAREFUL WHEN HANDLING AND PREPARING TUBES TO PREVENT BREAKAGE AND POSSIBLE INJURY. INSPECT **QBC** ACCUTUBES BEFORE USE. DO NOT USE CRACKED OR SCRATCHED TUBES. IF A TUBE BREAKS, REFER TO THE **QBC Autoread Plus** SYSTEM OPERATOR'S MANUAL FOR GUIDANCE.

ACCUTUBE STORAGE

AccuTubes are shipped in a vial. Store AccuTubes in the vial with stopper tightly in place. Remove tubes as needed and replace the stopper. Keep the desiccant strip in the vial. Store vials in a dry cupboard or drawer, at 60° to 90° F (16° to 32° C). Do not store in bright light.

An expiration date is printed on the vial. Do not use vials of tubes that are past the expiration date. Once the vial has been opened, tubes can be used for up to 30 days, as long as the expiration date has not passed. Use the space on the vial label to record the date the vial was opened.

The foam plug must be removed from newly opened AccuTube vials. Use the forceps provided with the **QBC** system.

AccuTubes must be visually inspected for correct position of stopper. The stopper must be partially seated.

BLOOD COLLECTION AND STORAGE

Refer to the **QBC Autoread Plus** System Operator's Manual for precautions and procedures related to blood collection.

AccuTubes can be used with blood samples collected either by venipuncture (venous blood) or by skin puncture (capillary blood). Venipuncture is the preferred collection method for hematology testing.^{4,5,6} Follow accepted practices when taking the sample.⁵

PREPARATION AND TESTING

Users should review the **QBC Autoread Plus** System Operator's Manual before using these procedures.

- Materials Provided: Catalog No. 423406, 100 AccuTube tests - 10 x 10 AccuTube SoftPack, with floats and instructions for use.
- Materials Required: **QBC Autoread Plus** System, lint-free tissues.
- Optional Materials: **QBC** Pipetter with AccuTube spacer, AccuTube Barcode Labeler.

Sample Storage and Timing

Venous blood samples may be stored at room temperatures 68° to 77° F (20° to 25° C) for up to 8 hours prior to preparation of AccuTubes.

Samples that can not be tested immediately must be refrigerated if the room temperature is above 77°F (25° C). Refrigerated samples stored at 36° to 46° F (2° to 8°C) are stable for up to 8 hours. Bring samples back to room temperature before you prepare an AccuTube.

The AccuTube should be centrifuged within 15 minutes after the plastic float is inserted.

Once tubes are centrifuged, the AccuTube is stable for 4 hours before testing, provided the AccuTube is stored vertically with stopper end down, and away from heat and intense light.

Pre-Test Conditions

Tests must be run at room temperature [68° to 98° F (20° to 37° C)].

The **QBC Autoread Plus** analyzer must be turned on and have passed the daily check rod test. Refer to the Operator's Manual for analyzer start-up and quality control procedures.

RUNNING A PATIENT SAMPLE

Step 1: Fill the AccuTube

NOTE: Do not allow the blood to touch the AccuTube stopper while performing this step.

Venous Blood

Gently mix the sample at least 8 times by inversion, or for 5 minutes on a mechanical mixer immediately before filling the AccuTube.

Tilt the blood tube, and place the open end of the AccuTube in contact with the blood.

Fill the AccuTube until the blood level is between the two black lines.

Wipe the outside of the AccuTube with lint-free tissue.

Note: For instructions on use of the pipetter provided with the **QBC Autoread Plus** System, refer to the Operator's Manual.

Capillary Blood

Place the open end of the AccuTube in contact with the finger puncture blood.

Hold tube close to horizontal to avoid air bubbles.

Fill the AccuTube until the blood level is between the two black lines.

Wipe the outside of the AccuTube with lint-free tissue.

Step 2: Rock the AccuTube to Mix

NOTE: Do not allow the blood to touch the AccuTube stopper while performing this step.

If blood stops mixing, loosen or remove stopper and finish mixing. Re-insert stopper.

Hold the AccuTube in the center and allow the blood to move toward the stopper end. Rock the AccuTube at least 5 times to mix well with orange coating. For each complete rock, allow the blood to move from the stopper end to the open end and back to the stopper end.

Step 3: Seat the Stopper

Turn the AccuTube upright (stopper end down). At this point the blood can touch the stopper.

Hold the AccuTube near the stopper end.

Seat the stopper by pressing down firmly on a hard surface.

Step 4: Insert the Float

NOTE: Do not touch the floats with your fingers.

Slide the open end of the AccuTube over a float until the float is partially inserted. Gently lift the AccuTube up and out of the SoftGrip. Push float against the back of the package to fully insert.

Step 5: Apply Bar Code Labels (Optional)

NOTE: The bar code labels come in pairs. The small label goes on the AccuTube, and the larger label goes on the patient record.

Make sure the labels are legible before you apply them to the tube and patient record.

Apply Label to AccuTube

Place the AccuTube into the slot at the top front portion of the labeler, with the stopper to the left.

Advance the thumb wheel until a small label is fed up to the top of the AccuTube.

Using thumb and forefinger, pick up the label and AccuTube together. Wrap the label completely around the AccuTube avoiding wrinkles or bubbles.

Apply Label to Patient Record

Advance the thumb wheel until the large label is exposed. Remove the label and place on the patient record.

Step 6: Centrifuge the AccuTube

Open the centrifuge lid and remove the metal cover.

Hold the AccuTube so the stopper is toward the outside edge of the rotor. Place the AccuTube in the centrifuge.

Balance the rotor by placing an empty AccuTube in the slot opposite the patient AccuTube. Place the removable metal cover over the rotor and screw in place until finger tight.

Close the lid and press ON/OFF button. The centrifuge will spin for 5 minutes.

The lid will unlatch when the centrifuge stops. Open the lid, unscrew the metal cover, and remove the AccuTubes.

Step 7: Analyze the Sample

Select the CBC mode on the **QBC Autoread Plus** Analyzer.

Open the analyzer door and place the AccuTube in the analyzer with the stopper to the left.

Close the analyzer door to start the test.

Step 8: Obtain Results

When the test is complete, the results will be displayed on the analyzer and printed.

Test results are cleared from the display when the analyzer door is opened. Do not open the analyzer door until the results are reviewed and additional printouts made.

Open the analyzer door and remove the sample.

Discard the sample in a biohazardous sharps container.

TECHNICAL INFORMATION

See the **QBC Autoread Plus** System Operator's Manual for information concerning accuracy, precision and operating range of **QBC** AccuTube hematology tests.

Interfering Substances

Hemolysis: Do not perform tests on visibly hemolyzed blood specimens.

Bilirubin: No effects on test results have been observed at bilirubin concentrations up to 20 mg/dL.⁷

Triglycerides: No effects on test results have been observed at triglycerides concentrations up to 1,800 mg/dL.⁷

Coumadin: Anticoagulant therapy has been shown to have no clinically significant effect on performance.⁷

Doxorubicin: Treatment with the anthracyclenic drug Doxorubicin does not appear to interfere with the **QBC** test method.⁷

Other drugs: The effects of other potentially interfering drugs and their metabolites on **QBC** tests have not been established.^{8,9}

QUALITY CONTROL

Instrument performance is checked daily by running the check rod. See the **QBC Autoread Plus** System Operator's Manual for analyzer start-up and quality control procedures.

LIMITATIONS

Quality medical care requires that laboratory values be correlated with each patient's symptoms and signs by a trained practitioner.

The Performance section of the **QBC Autoread Plus** System Operator's Manual lists the validated upper and lower limits of the operating range. Values above and below these validated ranges should be confirmed by an alternate method.

The AccuTube has been formulated to provide optimum packing and layering of normal cells. In a small number of patients, however, the system cannot read certain parameters, and will not report a value. User errors in processing or use of outdated or inappropriately stored tubes can also result in non-reported results. Practitioners must not assume that unreported values are normal; further testing with an alternative method is essential.

Automated granulocyte and lymphocyte/monocyte differential counts cannot replace the conventional manual differential. Due to the grouping by density of the cell populations by the **QBC** test method, the system cannot discriminate between normal and abnormal cell types in disease states characterized by the presence of abnormal white cell types or nucleated red blood cells. If abnormal cell populations are suspected, verification of **QBC** test results or testing and diagnosis by alternative methods is essential.

The combined lymphocyte/monocyte count should not be used to test for lymphocytopenia in evaluating patients with known or possible immunodeficiencies. Further evaluation of lymphocyte/monocyte counts in relevant situations must include a manual differential and lymphocyte subset analysis.

The presence of abnormally sized platelets may lead to discrepancies between the **QBC** test method platelet count, which is based on platelet mass, and results obtained with an impedance counter, which are based on measurement of particle number.

BIBLIOGRAPHY

1. Wardlaw, S.C. and Levine, R.A.: "Quantitative Buffy Coat Analysis," *JAMA* 5:617-620 (1983).
2. Zucker, R.M. and Casse, B. (1966) "The Separation of Normal Human Leukocytes by Density and Classification by Size," *Blood*, 34:5, 591-600.
3. Jackson, J.F. (1961) "Supravital Blood Studies, Using Acridine Orange Fluorescence," *Blood*, 17:2643-649.
4. Koepke, John A., (1991) "Specimen Collection for Cellular Hematology," *In Practical Laboratory Hematology*, Koepke, John A. (Ed.), Churchill Livingstone, New York, pp. 33-35.
5. Nelson, Douglas A. and Morris, Michael W. (1984) "Basic Methodology," *In Clinical Diagnosis and Management by Laboratory Methods*, Henry, John B. (Ed.) 17 Edition, W.B. Saunders Company, Philadelphia, pp. 578-579.
6. National Committee for Clinical Laboratory Standards. Procedures for the Collection of Diagnostic Blood Specimens by Skin Puncture; Approved Standard (1991) NCLC document H4-A3 Villanova, Pennsylvania 19085.
7. Data on file at Becton Dickinson Primary Care Diagnostic, Sparks, MD 21152-0370.
8. Yount, D.S., Pestaner, L.C. and Gibberman, V. (1975) "Effects of Drugs on Clinical Laboratory Tests," *Clin. Chem.*, 21, 313D, 3454D, 346D, 390D, 391D, 392D.
9. Elking, M.P. and Kabat, H. (1968) "Drug Induced Modifications of Laboratory Test Values," *Am. J. Hosp. Pharm.*, 25, 485.

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A division of Woodley Equipment Company Ltd