# QBC® STAR Tube

L-001078 Rev. D

English:	pages	1– 4	Italiano:	pagine	11– 14
Français:	pages	5– 8	Español:	páginas	15– 18
Deutsch:	Seiten	8– 11			

U.S. Patent Numbers: 6,002,474; 6,074,883; 6,080,366.

Kontakt venligst den lokale QBC-repræsentant for at bestille instrukser på dansk.  
Contacte o seu representante local da QBC para obter as instruções de utilização em português.  
Kontakta närmaste QBC-representant för att erhålla instruktioner på svenska.

## INTENDED USE

STAR Tubes are used with the **QBC STAR™** analyzer to provide a diagnostic hematology profile. Venous or capillary blood can be tested. The following quantitative values are obtained from a single STAR tube: Hematocrit, hemoglobin, mean corpuscular hemoglobin concentration (MCHC), platelet count, white blood cell count, granulocyte count (% and number), and lymphocyte/monocyte count (% and number).

## PRINCIPLES OF THE PROCEDURE

When a tube of blood is placed in a centrifuge and spun, the different types of cells in the blood separate into layers from the heaviest to the lightest. A special coating in the STAR tube stains the different cell layers. These layers are made larger by a plastic float inside the tube. The **QBC STAR** analyzer measures the cell layers and converts the measurements into reported values.<sup>1,2,3</sup>

## SUMMARY AND EXPLANATION

The **QBC STAR** instrument is a small, portable, table top centrifugal hematology analyzer. The instrument uses specially designed tubes known as the STAR Tube.

Each tube has two black fill lines, a white line, reagents and a cap. The tube is filled by capillary action with 70 µL of either anticoagulated venous blood or capillary blood. The blood is mixed with the reagents by rocking the tube, the float is inserted as the tube is capped by placing the tube cap over the collection end of the tube. After seating the cap, the tube is placed in the analyzer for automatic processing. The **QBC STAR** displays and prints the results at the completion of the processing cycle.

## REAGENTS

Each tube contains a maximum of 6.5 µg Acridine Orange, 0.65 mg Potassium Oxalate, 1.5 USP units Sodium Heparin (anticoagulant), 0.33 mg Dipotassium EDTA (anticoagulant), agglutinating agents and stabilizers.

## TUBE STORAGE

Store in a dry cupboard or drawer, at 60 to 90°F (16 to 32°C). Do not store in bright light. The tubes are individually packaged in a foil-plastic blister that also contains a desiccant sachet. Do not open or remove the tube from the unitized package until needed.

An expiration date is printed on the box and the unitized package. Do not use tubes that are past the expiration date. The tube must be used after opening the unitized package; do not remove the tube from the unitized package until needed.

## BLOOD COLLECTION AND STORAGE

Refer to the **QBC STAR** System User's Manual for precautions and procedures related to blood collection. STAR Tubes can be used with blood samples collected either by venipuncture (venous blood) or skin puncture (capillary blood). Venipuncture is the preferred collection method for hematology testing.<sup>4,5,6</sup> Follow accepted practices when taking the sample.<sup>6</sup>

## WARNINGS

For *in vitro* Diagnostic Use.

**Acridine orange reagent may be toxic; do not ingest. Avoid contact with skin, eyes and clothing.**

**Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"<sup>7-10</sup> and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.**

**STAR Tubes consist of a glass tube that extends beyond the plastic sleeve prior to capping the tube. Be careful when handling and preparing tubes to prevent breakage and possible injury. Inspect STAR Tubes before use. Do not use cracked or broken tubes. If a tube breaks, refer to the QBC STAR System User's Manual for guidance.**

### **PREPARATION AND TESTING**

Users should review the **QBC STAR** System User's Manual before using these procedures.

- Materials Provided: Catalog No. 429625, 100 Tube tests or Catalog No. 429011, 20 Tube tests – individually packaged, and assembled with cap, floats and instructions for use.
- Materials Required But Not Provided: **QBC STAR** Instrument.

### **SAMPLE STORAGE AND TIMING**

Venous blood samples may be stored at room temperatures 68 to 77°F (20 to 25°C) for up to 8 h prior to preparation of STAR Tubes.

Samples that cannot be tested immediately must be refrigerated if the room temperature is above 77°F (25°C). Refrigerated samples stored at 36 to 46°F (2 to 8°C) are stable for up to 8 h. Bring sample back to room temperature before you prepare a STAR Tube.

The tube should be placed in the **QBC STAR** analyzer promptly after filling with blood and capping, and no later than 15 min after filling with blood.

- Pre-Test Conditions

Tests must be run at room temperature [68 to 98°F (20 to 37°C)]

The **QBC STAR** analyzer must be turned on and have passed the system self test before processing samples. Refer to the Operator's Manual for analyzer start-up and quality control procedures.

### **RUNNING A PATIENT SAMPLE**

#### **STEP 1: FILL THE TUBE**

NOTE: Refer to **QBC STAR** Quick Reference Guide for illustrations of these steps.

#### **ANTICOAGULATED VENOUS BLOOD**

Gently mix the sample at least 12 – 15 times by inversion before filling the tube.

Tilt the sample vial, and place the collection tip of the STAR Tube in contact with the blood tube. Do not allow the plastic sleeve to contact the blood.

Fill the tube beyond the first black line and up to the second black line.

#### **CAPILLARY BLOOD**

Place the collection tip of the STAR Tube in contact with the blood droplet at the finger puncture site. Do not allow the plastic sleeve to contact the blood.

Fill the tube beyond the first black line and up to the second black line.

#### **ROCK THE TUBE TO MIX**

**NOTE:** Do not allow the blood to touch the tube plug while performing this step.

If blood stops mixing, discard tube and fill a new tube.

Rock the tube back and forth at least 4 times to mix the blood with the orange coating. For each complete rock, allow the blood to move from the open end to the bottom (plug) end and back to the open end.

Hold the tube at an angle allowing the blood to move down the tube towards the center of the tube.

#### **STEP 2: CAP THE TUBE**

Remove the cap from the tube by pulling the cap straight off the bottom of the tube.

Place the cap over the collection end of the tube by guiding the glass end of the tube into the center of the cap.

Push the cap on firmly.

#### **STEP 3: READ THE TUBE**

Promptly insert the tube into the **QBC STAR** analyzer oriented with the top of the tube (capped end) pointing into the analyzer.

Close the analyzer door.

Press the “start” button.

When the test is complete, the results are automatically displayed and printed.

**DO NOT OPEN THE ANALYZER DOOR UNTIL THE RESULTS ARE REVIEWED AND PRINTED.**

Test results are cleared from the display when the analyzer door is opened.

Open the analyzer door and remove the sample. Discard the STAR Tube in a biohazard container.

## **TECHNICAL INFORMATION**

See the **QBC STAR** System User’s Manual for information concerning accuracy, precision and operating range of **QBC STAR** hematology tests.

## **INTERFERING SUBSTANCES**

Hemolysis: Do not perform tests on visibly hemolyzed blood specimens.

Bilirubin: No effects on test results have been observed at bilirubin concentrations up to 20 mg/dL.<sup>11</sup>

Triglycerides: No effects on test results have been observed at triglycerides concentrations up to 1,800 mg/dL.<sup>11</sup>

Coumadin: Anticoagulant therapy has been shown to have no clinically significant effect on performance.<sup>11</sup>

Doxorubicin: Treatment with the anthracyclic drug Doxorubicin does not appear to interfere with the **QBC™** test method.<sup>11</sup>

Other drugs: The effects of other potentially interfering drugs and their metabolites on **QBC** tests have not been established.<sup>12,13</sup>

## **QUALITY CONTROL**

Instrument performance is checked automatically, during every sample run, by a series of electronic, optical and mechanical measurements. Blood sample results are only displayed if all electronic quality control measurements pass factory set limits. Successful quality checks are documented on the printout with sample results. See the **QBC STAR** Instrument Operator’s Manual for further information concerning analyzer start-up and quality control procedures.

## **LIMITATIONS**

Quality medical care requires that laboratory values be correlated with each patient’s symptoms and signs by a trained practitioner.

The Performance section of the **QBC STAR** System User’s Manual lists the validated upper and lower limits of the operating range. Values above and below these validated ranges are not printed or displayed by the analyzer and such samples must be tested by an alternate method.

The tube has been formulated to provide optimum packing and layering of normal cells. In small numbers of patients, however, the instrument cannot read certain parameters, and will not report a value. User errors in processing or use of outdated or inappropriately stored tubes can also result in non-reported results. Practitioners must not assume that unreported values are normal; further testing with an alternative method is essential.

Automated granulocyte and lymphocyte/monocyte differential counts cannot replace the conventional manual differential. Due to the grouping by density of the cell populations by the **QBC** test method, the instrument cannot discriminate between normal and abnormal cell types in disease states characterized by the presence of abnormal white cell types or nucleated red blood cells. If abnormal cell populations are suspected, verification of **QBC** test results or testing and diagnosis by alternative methods is essential.

The combined lymphocyte/monocyte count should not be used to test for lymphocytopenia in evaluating patients with known or possible immunodeficiencies. Further evaluation of lymphocyte/ monocyte counts in relevant situations must include a manual differential and lymphocyte subset analysis.

The presence of abnormally sized platelets may lead to discrepancies between the **QBC** test method platelet count, which is based on platelet mass, and results obtained with an impedance counter, which are based on measurement of particle number.

#### AVAILABILITY

Cat. No.	Description
429011	<b>QBC STAR™</b> 20 STAR Tubes
429625	<b>QBC STAR™</b> 100 STAR Tubes

#### REFERENCES

1. Wardlaw, S.C. and Levine, R.A.: Quantitative Buffy Coat Analysis: JAMA (1983); 5:617-620.
2. Zucker, R.M. and Casse, B. (1966) The Separation of Normal Human Leukocytes by Density and Classification by Size; Blood, 34:5, 591-600.
3. Jackson, J.F. (1961): Supravital Blood Studies, Using Acridine Orange Fluorescence; Blood, 17:2643-649.
4. Koepke, John A., (1991): Specimen Collection for Cellular Hematology; In Practical Laboratory Hematology, Koepke, John A. (Ed.) Churchill Livingstone, New York, pp. 33-35.
5. Nelson, Douglas A. and Morris, Michael W. (1984): Basic Methodology, In Clinical Diagnosis and Management by Laboratory Methods, Henry, John B. (Ed.) 17 Edition, W.B. Saunders Company, Philadelphia, pp. 578-579.
6. National Committee for Clinical Laboratory Standards. Procedures for the Collection of Diagnostic Blood Specimens by Skin Puncture; Approved Standard (1991) NCCLS document H4-A3 Villanova, Pennsylvania 19085.
7. National Committee for Clinical Laboratory Standards. 1997. Approved Guideline M29-A. Protection of laboratory workers from instrument biohazards and infectious disease transmitted by blood, body fluids and tissue. NCCLS, Wayne, Pa.
8. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80
9. U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
10. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
11. Data on file at QBC Diagnostics, Philipsburg, PA 16866.
12. Yount, D.S., Pestaner, L.C. and Gibberman, V. (1975): Effects of Drugs on Clinical Laboratory Tests; Clin, Chem., 21, 313D, 345D, 346D, 390D, 391D, 392D.
13. Elking, M.P. and Kabat, H. (1968) "Drug Induced Modifications of Laboratory Test Values," Am. J. Hosp. Pharm., 25, 485.

QBC, QBC Logo and QBC Diagnostics are trademarks of QBC Diagnostics, Inc. © 2004 QBC Diagnostics, Inc.

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Port Matilda, PA 16870 U.S.A.
A division of Woodley Equipment Company Ltd