

BOOST 4+ Flex **LARGE CAPACITY CENTRIFUGE**

Part # 00-389-009-000

Enjoy the highest performance and capacity with the BOOST 4+ Flex. Meet demanding processing needs with speeds up to 4000 RPM. Easily accommodate large tubes in four swing-out buckets.

Fully Customizable

The BOOST 4+ Flex offers a variety of extensive customization options. Personalize the front panel label and LED lid lights to match your unique branding, or pre-program custom RPM and g-force settings to meet your preferred processing specifications.

Spins Your Cycles

Program and recall up to 10 preset cycles with the BOOST 4+ Flex's backlit digital display. Can be set to run preset cycles only, or may be locked to just one preset cycle if desired.

Fits Your Tubes

Experience unmatched convenience with the BOOST 4+ Flex's unique insert system. Simply select the correct adapter for each tube size and switch them out as needed to spin every tube.

Safe and Durable

Enjoy years of worry-free operation with BOOST's imbalance protection and shatterproof lid, which can only be opened when the rotor is stopped - even if the unit loses power. Rigorous electrical testing ensures our centrifuges meet hospital equipment safety standards.


A Great Value

Unbeatable quality and flexibility for any budget. Each BOOST 4+ Flex centrifuge is ready to run with swing-out horizontal rotor, buckets, and 2-year warranty included in the affordable price.



Unique Rotor Design

Our patented high-volume horizontal rotor and custom inserts designed in-house ensure a perfect fit for any tube.



Applications	Max Speed	Max Capacity	Included	Warranty
Variable Tube Sizes Custom Disposables	4000 RPM 2950 xg RCF	4 x up to 1.9 x 5.4 in (4.8 x 13.7 cm) or 140 grams	Horizontal Rotor 2 x BOOST Bucket for Molded Rotor	2 years Extended warranty options available

Designed, built, and supported in the USA

BOOST 4+ Flex



Free 30-Day Evaluation (U.S. Only)

Part Numbers

		Part Numbers
	Centrifuge	00-389-009-000
Included	BOOST Bucket for Molded Rotor (2 qty)	03-1-0007-0123K
	.68" BOOST Bucket Adapter, Red (2 qty)*	03-1-0007-0112K
Optional	.85" BOOST Bucket Adapter, Red (2 qty)*	03-1-0007-0113K
	1.03" BOOST Bucket Adapter, Red (2 qty)*	03-1-0007-0114K
	1.2" BOOST Bucket Adapter, Red (2 qty)*	03-1-0007-0115K
	1.34" BOOST Bucket Adapter, Red (2 qty)*	03-1-0007-0116K
	1.56" BOOST Bucket Adapter, Red (2 qty)*	03-1-0007-0117K
	1 yr Warranty Extension	03-0-0007-0003
	2 yr Warranty Extension	03-0-0007-0004
3 yr Warranty Extension	03-0-0007-0005	

*also available in black adapters

Specifications

Sound Level	75 dBA
Max Speed RPM	4000 RPM
Max G-Force	2950 xg
Dimensions (W x D x H)	15 x 17 x 9 in 38.1 x 43.2 x 22.9 cm
Weight	39 lbs (17.69 kg)
Power	100-240 V – 50/60 Hz
Resting Tube Angle	90° (vertical)
Tube Spin Angle	90° (horizontal)
Maximum Load	140 g per bucket



druckerdiagnostics.com
sales@druckerdiagnostics.com

+1-866-265-1486 (U.S. only)
+1-814-692-7661