



DASH Centrifuges

High Performance STAT Processing

| Model | Rotor Type | Interface | Preset Cycles Time / Speed / RCF | Max Speed RPM / RCF | Up to 13 x 75 mm or 5 mL | Dimensions (W x D x H) |
|-------------------------------|-------------|-----------|---|------------------------|--------------------------|---------------------------------|
| DASH Apex 4 00-053-009-000 | Fixed Angle | 3-Setting | 30 sec / 8500 RPM / 4440 xg 2 min / 8500 RPM / 4440 xg 3 min / 8500 RPM / 4440 xg | 8500 / 4440 | 4 | 7 x 7 x 7 in 18 x 18 x 18 cm |

| Model | Rotor Type | Interface | Preset Cycles Time / Speed / RCF | Max Speed RPM / RCF | Up to 16 x 100 mm or 10 mL | Up to 16 x 125 mm or 15 mL | 50 mL | Dimensions (W x D x H) |
|--------------------------------|-------------|------------------------|---|------------------------|----------------------------------|----------------------------------|-------|-----------------------------------|
| DASH Coag 00-676-009-001 | Fixed Angle | 3-Setting | 2 min / 6600 RPM / 4400 xg 3 min / 6600 RPM / 4400 xg 5 min / 5200 RPM / 2700 xg | 6600 / 4400 | 12 | - | - | 12 x 14 x 9 in 30 x 36 x 23 cm |
| DASH Apex 6 00-076-009-004 | Horizontal | 3-Setting | 3 min / 5300 RPM / 4000 xg 5 min / 5300 RPM / 4000 xg 7 min / 3800 RPM / 2000 xg | 5300 / 4000 | 6 | - | - | 12 x 14 x 9 in 30 x 36 x 23 cm |
| DASH Flex 6 00-176-009-000 | Horizontal | 10-Setting, Digital | 3 min / 5300 RPM / 4000 xg 5 min / 5300 RPM / 4000 xg 10 min / 3800 RPM / 2000 xg | 5400 / 4150 | 6 | - | - | 12 x 14 x 9 in 30 x 36 x 23 cm |
| DASH Apex 12 00-083-009-003 | Horizontal | 3-Setting | 3 min / 5200 RPM / 4000 xg 5 min / 5200 RPM / 4000 xg 7 min / 3600 RPM / 2000 xg | 5200 / 4000 | 12 | - | - | 13 x 15 x 9 in 33 x 38 x 23 cm |
| DASH Flex 12 00-183-009-000 | Horizontal | 10-Setting, Digital | 3 min / 5200 RPM / 4000 xg 5 min / 5200 RPM / 4000 xg 7 min / 3600 RPM / 2000 xg | 5200 / 4000 | 12 | - | - | 13 x 15 x 9 in 33 x 38 x 23 cm |
| DASH Apex 24 00-084-009-005 | Horizontal | 3-Setting | 5 min / 4300 RPM / 3150 xg 6 min / 4300 RPM / 3150 xg 7 min / 3500 RPM / 2100 xg | 4300 / 3150 | 24 | 12 | 6 | 15 x 17 x 9 in 38 x 43 x 23 cm |
| DASH Flex 24 00-184-009-000 | Horizontal | 10-Setting, Digital | 5 min / 4300 RPM / 3150 xg 6 min / 4300 RPM / 3150 xg 7 min / 3500 RPM / 2100 xg | 4500 / 3450 | 24 | 12 | 6 | 15 x 17 x 9 in 38 x 43 x 23 cm |

(adapters available for a variety of tube sizes)



SERO 12 - Blood Banking Centrifuge

A drop-in replacement for the BD SeroFuge

| Model | Rotor Type | Interface | Preset Cycle Time / Speed / RCF | Max Speed RPM / RCF | 12 Place Rotor Up to 75mm | 6 Place Rotor Up to 100 mm | Dimensions (W x D x H) |
|---------------------------|-------------|------------------------|------------------------------------|------------------------|------------------------------|-------------------------------|-----------------------------------|
| SERO 12 00-483-009-001 | Fixed Angle | 10-Setting, Digital | 15 sec / 3600 RPM / 1040 xg | 3600 / 1040 | 6 capped 12 uncapped | 6 | 13 x 15 x 9 in 33 x 38 x 23 cm |



All centrifuges are available
for a free 30-day evaluation.
Visit: druckerdiagnostics.com/evaluation-program


HORIZON Centrifuges
Quality Routine Chemistry

| Model | Rotor Type | Interface | Preset Cycles Time / Speed / RCF | Max Speed RPM / RCF | Up to 16 x 100 mm or 10 mL | Up to 16 x 125 mm or 15 mL | 50 mL | Dimensions (W x D x H) |
|--|-------------|------------------------|--|------------------------|----------------------------------|----------------------------------|-------|-----------------------------------|
| HORIZON 6 00-276-009-000 | Horizontal | 3-Setting | 10 min / 3600 RPM / 1800 xg 15 min / 3300 RPM / 1500 xg 5 min / 1900 RPM / 500 xg | 3800 / 2000 | 6 | - | - | 12 x 14 x 9 in 30 x 36 x 23 cm |
| HORIZON 6 FA 00-276-009-001 | Fixed Angle | 3-Setting | 10 min / 3900 RPM / 1600-1850 xg 15 min / 3800 RPM / 1500 xg 5 min / 2200 RPM / 500-600 xg | 3900 / 1850 | 6 | 6 | - | 12 x 14 x 9 in 30 x 36 x 23 cm |
| HORIZON 6 Flex 00-376-009-000 | Horizontal | 10-Setting, Digital | 10 min / 3600 RPM / 1800 xg 15 min / 3300 RPM / 1500 xg 5 min / 1900 RPM / 500 xg | 3800 / 2000 | 6 | - | - | 12 x 14 x 9 in 30 x 36 x 23 cm |
| HORIZON 6 Flex FA 00-376-009-001 | Fixed Angle | 10-Setting, Digital | 10 min / 3900 RPM / 1600-1850 xg 15 min / 3800 RPM / 1500 xg 5 min / 2200 RPM / 500-600 xg | 3900 / 1850 | 6 | 6 | - | 12 x 14 x 9 in 30 x 36 x 23 cm |
| HORIZON 12 00-238-009-000 | Horizontal | 3-Setting | 10 min / 3500 RPM / 1800 xg 15 min / 3200 RPM / 1500 xg 5 min / 1800 RPM / 500 xg | 3700 / 2000 | 12 | - | - | 13 x 15 x 9 in 33 x 38 x 23 cm |
| HORIZON 12 Flex 00-383-009-000 | Horizontal | 10-Setting, Digital | 10 min / 3500 RPM / 1800 xg 15 min / 3200 RPM / 1500 xg 5 min / 1800 RPM / 500 xg | 3700 / 2000 | 12 | - | - | 13 x 15 x 9 in 33 x 38 x 23 cm |
| HORIZON 24 00-284-009-000 | Horizontal | 3-Setting | 10 min / 3300 RPM / 1800 xg 15 min / 3000 RPM / 1500 xg 5 min / 1700 RPM / 500 xg | 3400 / 2000 | 24 | 12 | 6 | 15 x 17 x 9 in 38 x 43 x 23 cm |
| HORIZON 24 Flex 00-384-009-000 | Horizontal | 10-Setting, Digital | 10 min / 3300 RPM / 1800 xg 15 min / 3000 RPM / 1500 xg 5 min / 1700 RPM / 500 xg | 3400 / 2000 | 24 | 12 | 6 | 15 x 17 x 9 in 38 x 43 x 23 cm |

(adapters available for a variety of tube sizes)


600 Series Centrifuges
Compact Routine Chemistry

| Model | Rotor Type | Interface | Preset Cycle Time / Speed / RCF | Max Speed RPM / RCF | Up to 16 x 100 mm or 10 mL | Up to 16 x 125 mm or 15 mL | 50 mL | Dimensions (W x D x H) |
|------------------------|-------------|-----------------|------------------------------------|------------------------|----------------------------------|----------------------------------|-------|-----------------------------------|
| 614B 00-004-009-001 | Fixed Angle | 1-Setting, Dial | 1-30 min / 3250 RPM / 1200 xg | 3250 / 1200 | 6 | 6 | - | 12 x 14 x 9 in 30 x 36 x 23 cm |
| 642B 00-016-009-001 | Horizontal | 1-Setting, Dial | 1-30 min / 3380 RPM / 1600 xg | 3380 / 1600 | 6 | - | - | 12 x 14 x 9 in 30 x 36 x 23 cm |
| 642E 00-018-009-001 | Horizontal | 1-Setting | 10 min / 3380 RPM / 1600 xg | 3380 / 1600 | 6 | - | - | 12 x 14 x 9 in 30 x 36 x 23 cm |
| 642M 00-077-009-000 | Horizontal | 1-Setting | 10 min / 3400 RPM / 1600 xg | 3400 / 1650 | 6 | - | - | 12 x 14 x 9 in 30 x 36 x 23 cm |

(adapters available for a variety of tube sizes)

Questions? Contact Drucker Diagnostics Now.