

Announcing the TrueBond Duet, our latest and most compact centrifuge for animal health professionals.



Why choose TrueBond Duet instead of StatSpin™ VT?

- **Excellent Quality & Value.** TrueBond Duet combines durable construction with affordability, delivering reliable performance at a great price.
- **Compact & Capable.** TrueBond Duet is our most compact centrifuge for animal health professionals and comes with two rotors to support a variety of tube sizes.
- **Warranty.** TrueBond Duet carries a 2-year overall warranty and calibration guarantee. We also offer a lifetime warranty on all rotors.
- **Customer Service & Tech Support.** TrueBond Duet centrifuges are designed, built, and supported in the USA. Our customer service team is always just a phone call away.

The TrueBond Duet includes a 2-place, fixed angle rotor and a microhematocrit rotor to support the centrifugation of multiple tube sizes.



TrueBond Duet



StatSpin™ VT

Interface	5-Setting	5-Setting
Preset Cycles Time / RPM / RCF	30 sec / 15800 RPM / 12000 xg 1.5 min / 15800 RPM / 12000 xg 2.5 min / 15800 RPM / 12000 xg 45 sec / 9800 RPM / 3900 xg 2 min / 15800 RPM / 13700 xg	30 sec / 15800 RPM / 12000 xg 90 sec / 15800 RPM / 12000 xg 150 sec / 15800 RPM / 12000 xg 45 sec / 9800 RPM / 3900 xg 120 sec / 15800 RPM / 13700 xg
Max Speed RPM / RCF	15800 / 13700 (Microhematocrit) 15800 / 12000 (1.5 - 2 mL)	15800 / 13700
1.5 - 2 mL	2 (Fixed Angle Rotor)	2 (RT12 Rotor)
Microhematocrit (Up to 40 mm)	12 (MHCT Rotor)	12 (RH12 Rotor)
Dimensions (W x H x D)	7 x 7 x 7 in 18 x 18 x 18 cm	6 x 6.6 x 8.6 in 15.2 x 16.8 x 21.8 cm

A centrifuge designed for animal health professionals




Why choose TrueBond instead of LW E8?

- **Easy to Use & Monitor.** TrueBond comes equipped with a fully digital, user-friendly interface and LED lid lighting to indicate cycle status and is designed to spin blood, urine, fecal and microhematocrit tubes.
- **Excellent Quality & Value.** TrueBond is price-favorable to the lab and offers a model for every budget.
- **Warranty.** TrueBond carries a 2-year overall warranty and calibration guarantee. We also offer a lifetime warranty on all rotors.
- **Customer Service & Tech Support.** TrueBond centrifuges are designed, built, and supported in the USA. Our customer service team is always just a phone call away.



TrueBond



LW E8 Touch

Interface	10-Setting	10-Setting
Preset Cycles Time / RPM / RCF	10 min / 3900 RPM / 1600 xg 5 min / 2000 RPM / 500 xg 5 min / 1600 RPM / 300 xg 6 min / 3800 RPM / 1500 xg	10 min / 3500 RPM / 1534 xg 5 min / 1800 RPM / 406 xg 6 min / 1500 RPM / 282 xg 6 min / 3500 RPM / 1534 xg
Max Speed RPM / RCF	3900 / 1850	3500 / 1534
Microhematocrit	8	8
Up to 100 mm or 10 mL	6	8
Up to 125 mm or 15 mL	6	8
Dimensions (W x H x D)	12 x 14 x 9 in 30 x 36 x 23 cm	11 x 9.5 x 11 in 28 x 24 x 28 cm



All centrifuges are available for a free 30-day evaluation.
Learn more at: druckerdiagnostics.com/evaluation-program