



TrueBond Set 3

Set-and-Lock Veterinary Centrifuge

Part # 00-276-009-007

Using the TrueBond Set 3 set-and-lock veterinary centrifuge, vet laboratories can now spin microhematocrit tubes, fecal, blood, urine, and more in one centrifuge.

Spin Any Tube

TrueBond Set 3 is designed to spin all your tubes, from microhematocrit to urine and everything in between. Includes all necessary adapters.

Easy to Use

Minimize training time with the user-friendly interface. The Truebond Set 3 features simple controls designed for maximum repeatability and three presets for all common veterinary needs.

Easy to Monitor

Process samples more efficiently with TrueBond Set 3's LED lid lighting. Once the unit is started, the display shows remaining spin

Safe and Durable

Enjoy years of safe, worry-free operation with the brushless motor, imbalance protection, and a clear, shatterproof lid, which locks until the centrifuge has stopped completely.

A Great Value

Each TrueBond Set 3 centrifuge is ready to run with a fixed angle rotor, two sets of tube carriers, microtainer adapters, microhematocrit adapters and reader card, and a 2-year warranty included in the competitive price.



Applications	Max Speed	Max Capacity	Included	Warranty
Blood, Urine, Fecal, Microhematocrit	3900 Max RPM 1850 xg Max RCF	6x blood, urine, and fecal tubes Up to 16 x 125 mm or 15 mL 8x microhematocrits	Rotor (fixed angle) 6x 3-10 mL blue tube carriers 6x 15 mL orange tube carriers 6x .5-1 mL adapters Microhematocrit adapters and reader card	2 years Rotor: lifetime warranty Extended warranty options available

Free 30-Day Evaluation (U.S. Only)



Use the three presets for common veterinary applications, or program your validated cycles

Simple two-button controls

Lid lights indicate cycle status: ready, running, done

QR code for easy access to manual and more



Standard Preset Cycles		Blue Tube Carriers (Up to 16 x 100mm/10 mL tubes)	Orange Tube Carriers (Up to 16 x 125 mm/15 mL tubes)
Cycle 1	10 min @ 3900 RPM	1650 xg	1850 xg
Cycle 2	5 min @ 2000 RPM	425 xg	500 xg
Cycle 3	6 min @ 3800 RPM	1550 xg	N/A

Part Numbers	
Centrifuge	00-276-009-007
75-100 mm Blue Carriers (6 qty)	7713079K
125 mm Orange Carriers (6 qty)	7713044K
0.5-1 mL Purple Adapters (6 qty)	03-1-0007-0059K
Microhematocrit Adapters (2 qty)	02-004-0-0027K
Veterinary Microhematocrit Reader Card	02-0006-0-0046
1 Year Warranty Extension	03-0-0007-0001
2 Year Warranty Extension	03-0-0007-0002
3 Year Warranty Extension	03-0-0007-0006

Specifications	
Sound Level	64 dBA
Max RPM	3900 RPM
Max G-Force	1850 xg
Dimensions (W x D x H)	12 x 14 x 9 in 30 x 36 x 23 cm
Weight	12 lbs (5.4 kg)
Power	100 to 240 V – 50/60 Hz



druckerdiagnostics.com
sales@druckerdiagnostics.com

+1-866-265-1486 (U.S. only)
+1-814-692-7661

Protected by U.S. Patents #6,811,531. #D718,463S #D734,489S Other Patents Pending.