

Complete Fluorescence Systems

424350

QBC Malaria Test System with ParaLens Advance™, Microscope, and Centrifuge

The QBC Malaria Test System with ParaLens Advance, Microscope, and Centrifuge is a complete transmission light microscope and epi-fluorescence microscopy system for malaria parasite staining and review. Includes research grade microscope with 4x, 10x, 40x and 100x objectives, ParaLens Advance microscope attachment with 60x objective, QBC Capillary Centrifuge (120/230 V, 50/60 Hz operation), and the QBC Malaria Test (Box of 2000).

424450

QBC Malaria Test System with ParaLens Advance and Centrifuge

The QBC Malaria Test System with ParaLens Advance and Centrifuge is an epi-fluorescence microscopy system for malaria parasite staining and review. Includes ParaLens Advance microscope attachment with 60x objective, QBC Capillary Centrifuge (120/230 V, 50/60 Hz operation), and the QBC Malaria Test (Box of 2000).

Test Kits

253037

QBC Malaria Test (Box of 100)

The QBC Malaria test is a capillary tube coated with anticoagulants and acridine orange used to stain malaria parasites for detection using fluorescence microscopy. Includes 10 individual trays of 10 tests, with precision plastic floats and closures.

253005

QBC Malaria Test (Box of 2000)

Includes 8 individual boxes of 250 tests, with precision plastic floats and closures.

Accessories and Consumables

424278

ParaViewer

The ParaViewer is a plastic capillary tube holder, used to review the QBC Malaria Test on a microscope stage.

424226

QBC Workstation

Sturdy plastic workstation that can be used to hold 100 tubes of varying sizes, including QBC Malaria Tubes.

427708

Immersion Oil (Type B)

This 30 mL bottle of immersion oil is needed to view the QBC Malaria Test under 60x magnification.

427752

Immersion Oil (Type B) (Pack of 4)

420000

1.8 mm Blade Lancets (Box of 100)

These spring loaded lancets are ideal for drawing capillary blood for use with QBC Tubes. The 1.8 mm version is designed for pediatric use.

420001

2.3 mm Blade Lancets (Box of 100)

Same as p/n 420000 except for use on adult patients.

Related Equipment

For a complete listing of related accessories for the ParaLens Advance or the QBC Capillary Centrifuge, please refer to their respective Product Description Lists.

