


HORIZON

HORIZON 6, 12, & 24 Flex PROGRAMMABLE *ROUTINE* CENTRIFUGES



Standardize the whole lab on these centrifuges to simplify work flows, training, warranty, and service. Move them around to address bottle necks.

Tubes	Applications	Centrifugation	Included	Warranty
6, 12, or 24 75 & 100 mm	Plasma Serum Urine	Horizontal Max 2,000 xg (RCF)	Rotor Buckets or tube holders	2 years Rotor: lifetime warranty Extended warranty available

Designed, built, and supported in the USA



Horizon 6 FLEX



Horizon 12 FLEX




Horizon 24 FLEX

Part Numbers	Horizon 6	Horizon 12	Horizon 24
Centrifuge	00-276-009-000	00-283-009-000	00-284-009-000
Replacement 75 – 100 mm tube holder/Bucket	7713079	7713079	7713023
50 mL tube holder			7713017
Cap			7713035
SmartView Platform	00-079-009-001		
1 yr Extended Warranty	03-0-0007-0001	03-0-0007-0001	03-0-0007-0001
2 yr Extended Warranty	03-0-0007-0002	03-0-0007-0002	03-0-0007-0002

- Quiet operation
- Positive safety latch
- Maintenance free brushless motor
- Safety lid switch prevents centrifuge from running when lid is open
- Clear shatterproof lid

Cool Operation

Unique design of the motor housing and rotation chamber keeps specimens cool as required by CLSI H21-A5, even under the heaviest workloads.



Specifications	Horizon 6	Horizon 12	Horizon 24
Sound level	63.6 dB	55.4 dB	58.2 dB
Dimensions	8 x 11 x 13.5 in 20 x 28 x 34 cm	8.5 x 12.5 x 15.5 in 23 x 32 x 39 cm	9 x 14.5 x 17.0 in 23 x 37 x 43 cm
Weight	12 lbs (5.4. kg)	30 lbs (13.7. kg)	39 lbs (17.7. kg)
Power	100 to 240 V – 50 & 60 Hz		