

POINT-OF-CARE USE OF LED FLUORESCENCE MICROSCOPY COMBINED WITH ULTRASOUND IN THE DIAGNOSIS OF EXTRA-PULMONARY TB: PRELIMINARY RESULTS

Maria T. Giordani¹, Enrico Brunetti², Michela Pascarella³, Tom Heller⁴, Walter F. Kuhn⁵, ¹Infectious and Tropical Diseases Unit, San Bortolo Hospital, Vicenza, Italy, ²Division of Infectious and Tropical Diseases, University of Pavia, IRCCS San Matteo Hospital Foundation, Pavia, Italy, ³Microbiology Department, San Bortolo Hospital, Vicenza, Italy, ⁴Klinikum Perlach, Internal Medicine, Muenchen, Germany, ⁵Center for Operational Medicine, Georgia Health Sciences University, Augusta, GA, United States

INTRODUCTION

Extrapulmonary TB (EPTB) is difficult to diagnose and early diagnosis and therapy mean better prognosis and higher survival rates particularly in patients co-infected with HIV.¹ Ultrasound (US) guided fine needle aspiration is crucial in abdominal EPTB in order to obtain samples for direct microscopy and culture, which remains the reference standard.² LED Fluorescence Microscopy (LED-FM) is a potential alternative to conventional fluoroscopy in the identification of Acid Fast Bacilli (AFB).³

PATIENTS AND METHODS

- 10 patients, 5 m/ 5 f, age 26-72 y.o. (median 33.5±12.3), 4 of them HIV+, 2 with diabetes, 1 on peritoneal dialysis.
- All patients underwent US-guided aspiration with Chiba needle (*BIOPSYBELL SRL Mirandola -MO- Italy*) to obtain sample for suspected TB
- Tests done on the pus: LED-fluorescence microscopy (LED-FM) (*QBC diagnostic, Port Matilda, PA USA*), Auramine-Rhodamine fluorescence microscopy (AR-FM), Zhiel-Nielsen (Z-N), PCR (*DNA SDA BD Probe Tec, Franklin Lakes, NJ USA*), culture (broth)
- 10 samples of culture-proven negative pus were used as control group

AIM OF THE WORK

To compare traditional diagnostic tests used in our laboratory for TB diagnosis and LED-FM in patients with EPTB, presenting pus collections available for US-guided aspiration

RESULTS




Abdominal abscess at the mesenteric root



Paravertebral abscess in Pott's disease




Tubercular osteomyelitis



LED-FM



Different concentrations of mycobacteria in parafilm microscopic field (40 x)



AR stain and Z-N are the standard tests in our Lab

| PATIENTS | Examined sample | LED-FM | AR | Broth culture/PCR | Final diagnosis |
|--|---|---|---|--|--|
| A.M, f, 30 y.o., Nigeria, HIV1 + | Pus from US guided aspiration of a giant neck mass | ++++ | ++++ | MDR-MBT/MBT | EPTB due to disseminated MDR-TB in AIDS, CRF |
| F.C.,m, 43 y.o., Ivory Coast, HIV1 + | US guided FNAB from the abdominal lymphnodes (figure) | +, many artifacts | Negative China ink positive for Histoplasmosis | Neg/Neg Culture positive for Histoplasmosis | Histoplasmosis in AIDS with possible TB reactivation |
| P. B.,m, 32 y.o., Ghana, HIV1 and 2 + | Pus (FNAB) from axillary lymphnodes | +++ | +++ | MBT/MBT | Disseminated TB in AIDS |
| Y.B.,f, 33 y.o., Burkina Faso, HIV1 + | Pus from US guided aspiration from abd LN | ++++ | ++++ | MBT/MBT | Disseminated TB in AIDS |
| J.M., m, 72 y.o., Italy | Pus from US-guided aspiration of the left LC abscess | ++ | Gram stain: G+ cocci. A-R negative. | MBT/MBT Culture pos. for CONS (<i>Staph. scheiffleri</i>) | Reactivated EPTB, bacterial superinfection, immunocompromised patient (diabetes, peritoneal diagnosis) |
| A.H.,f, 34 y.o., Morocco | US guided aspiration of pus from cervical lymphnodes | + | Negative. | MBT and MAC/not done | Double mycobacterial infection in immunocompetent patient. |
| R.N., f., 31 y.o., India | US guided drainage with pig tail cath 18 G (figure) | ++ | ++ | MBT/MBT | Bilateral paravertebral abscess in Pott's disease |
| B.T., m, 34 y.o., Bangladesh | US guided aspiration, surgical toilette of a giant mass in the neck | +++ | + / +++ (different samples) | MDR-MBT/MBT | Disseminated MDR-TB in decompensated diabetes |
| K.B., m, 35 y.o., India | US guided aspirated pus from the neck mass | + | + | MBT/MBT | Reactivated latent EPTB |
| E. K, f, 26 y.o., Tunisia | US-guided drainage with a pig tail cath 18 G | ++ | negative | MBT/MBT | Tubercular osteomyelitis (figure) |
| Total: 10 patients, 5 m/f 1/1, 26-72 y.o.(M 33.5±12.3) 4 HIV+ | | 10 positive (1 suspected false positive) | 7 positive/3 negative | 9 positive/1 negative for MBT | |

DISCUSSION AND CONCLUSIONS

LED-FM combined with US shows promise as a rapid point-of-care method for the diagnosis of EPTB, particularly in resource-limited settings where EPTB is highly endemic and laboratories are lacking^{4,5}. Studies on large series are needed to evaluate the sensitivity and specificity of this approach.

REFERENCES

1. Giordani MT, Brunetti E, Binazzi R, et al. Extrapulmonary mycobacterial infections in a cohort of HIV-positive patients: ultrasound experience from Vicenza, Italy. *Infection* 2012 Sep 24. [Epub ahead of print]
2. Heller T, Goblirsch S, Wallrauch C, Lessells R, Brunetti E, Abdominal tuberculosis: sonographic diagnosis and treatment response in HIV-positive adults in rural South Africa. *Int J Infect Dis.* 2010 Sep;14 Suppl 3:e108-12.
3. Kuhn W, Armstrong D, Atteberry S, Dewbrey E, Smith D, Hooper N. Usefulness of the parafilm fluorescent microscope adaptor for the identification of mycobacteria in both field and laboratory settings. *Open Microbiol J.* 2010 Apr 30;4:30-3
4. Alvarez-Uria G, Azcona JM, Midde M, Naik PK, Reddy S, Reddy R. Rapid Diagnosis of Pulmonary and Extrapulmonary Tuberculosis in HIV-Infected Patients. Comparison of LED Fluorescent Microscopy and the GeneXpert MTB/RIF Assay in a District Hospital in India. *Tuberculosis research and treatment* 2012;2012:932862.
5. Das D, Selvakumar N. Can LED fluorescence microscopy replace Ziehl-Neelsen microscopy in tuberculosis detection? [Correspondence]. *Int J Tuberc Lung Dis.* 2012 Nov;16(11):1558.