

Front-Load Lockbox



Can pay for itself in one day of bad weather

Heat and cold can both render specimens unusable. If a site collects six samples in a day and sets them outside for collection in a standard metal lockbox, extreme temperatures could ruin all six samples at a cost of up to \$150 in re-draws. With SampleSafe protecting the specimens, no re-draws are required - saving \$150, more than the cost of the SampleSafe, in one day.

Time to Overheat Comparison

100°F (38°C) Outside Temperature, with Cold Pack




Insulated Steel Box



SampleSafe™ Lockbox

Time to Freezing Comparison

15°F (-9°C) Outside Temperature



Insulated Steel Box



SampleSafe™ Lockbox

FRONT LOAD LOCKBOX

Height 14 in. / 35.56 cm
 Width 15 in. / 38.1 cm
 Depth 8 in. / 20.32 cm
 Weight 11.5 lbs. / 5.22 kg
 Internal Volume 570 in³ / 9,340.62 cm³
 Warranty 5 years




Front Load Lockbox
 P/N: 00-078-009-000

Also available:










Top Load Lockbox
 P/N: 00-071-009-000



Secure anti-theft mounting options (door hanging, under door, or security cable).

Lid can be customized with your logo.

OPTIONAL ACCESSORIES

ANTI-THEFT DEVICES		PART NUMBER
	Door Hanging Bracket - 24"	02-004-0-0007
	Door Hanging Bracket - 36"	00-078-109-000
	Under Door Bracket Kit	02-004-0-0002
	Anchor Bracket Kit (Use with Cable)	02-004-0-0003
	Anchor Cable (6 Foot)	03-1-0012-0054
	Wall Mounting Kit (Bracket & Screws)	02-004-0-0009
FOR ADDITIONAL PERFORMANCE		PART NUMBER
	Cold Pack (6" x 9" x .8")	03-1-0012-0053