

DASH APEX 6, 12, & 24 SET & LOCK STAT CENTRIFUGES

Discover a better way to spin STAT with Drucker DASH Apex STAT centrifuges. Minimize lab error, reduce turnaround time, and streamline training with the DASH Apex's simple controls and powerful performance.

Easy to Monitor

Prevent forgotten samples and reduce TAT with LED lid lighting (patent pending). Bright lights indicate when the centrifuge is ready to load (off), running (on), or done (flashing). Lights continue flashing until the lid is opened.

Easy to Use

The simple two-button interface makes it easy to select a preset or program your custom validated STAT cycles. Use all three presets, or set and lock to a single cycle for maximum repeatability. Spin 75-100 mm (3-10 mL) tubes with the same rotor and universal tube holders.

Eliminate Wait Time

Further shorten TAT by using a centrifuge array in accessioning to eliminate queuing. The array takes advantage of the DASH Apex's compact size to allow samples to be spun as soon as they arrive.

Safe and Durable

Enjoy years of worry-free operation with the DASH Apex's brushless motor and rugged, reinforced carbon fiber components. Enhance lab safety with imbalance protection and a clear, shatterproof lid, which locks until the centrifuge has stopped completely.

A Great Value

Each centrifuge comes ready to run with a horizontal rotor, universal tube holders, and 2-year warranty included in the affordable price.



Tube Capacity	6, 12, or 24 (75-100 mm/3-10 mL)
Applications	STAT chemistry, coag, PPP
Centrifugation	Horizontal (max 4,000 xg/5,500 RPM)
Included	Rotor, tube holders or buckets
Warranty *	2 years (rotor: lifetime warranty)

* *Extended warranty available*



DASH Apex 6



DASH Apex 12



DASH Apex 24

Designed, built, and supported in the USA

DASH APEX 6, 12, & 24



Free 30-Day Evaluation (U.S. Only)

Part Numbers	DASH Apex 6	DASH Apex 12	DASH Apex 24
Centrifuge	00-076-009-004	00-083-009-003	00-084-009-005
Replacement 75-100 mm tube holder/Bucket	03-1-0007-0046	03-1-0007-0046	03-1-0007-0046
Adapters, 0.5 to 1 mL	7713068 (1 qty)	7713068 (1 qty)	7713068 (1 qty)
Adapters, 1.5 to 2 mL	7713065 (1 qty)	7713065 (1 qty)	7713065 (1 qty)
1 yr Warranty Extension	03-0-0007-0003	03-0-0007-0003	03-0-0007-0003
2 yr Warranty Extension	03-0-0007-0004	03-0-0007-0004	03-0-0007-0004
3 yr Warranty Extension	03-0-0007-0005	03-0-0007-0005	03-0-0007-0005
* Custom Settings	00-100-100-101	00-100-100-101	00-100-100-101

* Add this part number to your centrifuge order to have your custom validated cycles pre-programmed at the factory

Factory Settings	DASH Apex 6			DASH Apex 12			DASH Apex 24		
	RPM	Time	G-Force	RPM	Time	G-Force	RPM	Time	G-Force
3 min (5 min for 24)	5,300	3	4,000	5,200	3	4,000	4,200	5	3,000
5 min (6 min for 24)	5,300	5	4,000	5,200	5	4,000	4,200	6	3,000
Custom	3,800	7	2,000	3,600	7	2,000	3,500	7	2,000

Specifications	DASH Apex 6	DASH Apex 12	DASH Apex 24
Sound level	61 dBA	61 dBA	64 dBA
Dimensions (W x D x H)	11 x 13.5 x 8 in 28 x 34 x 20 cm	12.5 x 14.75 x 9 in 32 x 37 x 23 cm	14.5 x 17 x 9 in 37 x 43 x 23 cm
Weight	12 lbs (5.4 kg)	34 lbs (15 kg)	39 lbs (17 kg)
Power	95 to 253 VAC – 50/60 Hz		

FDA LISTED



RoHS
Compliant



ISO13485
certified



druckerdiagnostics.com
sales@druckerdiagnostics.com

+1-866-265-1486 (U.S. only)
+1-814-692-7661

