

429001

QBC STAR™ Hematology Analyzer

The QBC STAR is a dry hematology analyzer that delivers a 9-parameter Complete Blood Count with the touch of a button. The STAR's unique microprocessor performs self-calibration, ensuring accuracy. The STAR comes with 2 rolls of thermal printer paper, screwdriver, forceps, patient ID keyboard, 10 each of 1.8 mm and 2.3 mm blade lancets, 8GB USB memory stick, and external power supply and cord. For 120 V, 60 Hz/230 V, 50 Hz operation.

Tubes

429625

QBC STAR Tubes (Box of 100)

A box of 100 STAR Tubes. STAR Tubes are the only tubes designed to work with the STAR.

429011

QBC STAR Tubes (Bag of 20)

Consumables & Accessories

424304

Hematology Controls

These control samples are designed to test blood cell count values, and are a necessity for U.S. customers to meet CLIA requirements for the STAR. The controls are shipped every three weeks, and contain 6 vials each of two different blood samples.

424305

Extended Range Controls

These control samples are designed to meet CLIA requirements for calibration verification. They contain 2 vials each of 3 different blood samples. Note: Unlike hematology controls (424304), extended range controls are available on an as-needed basis.

420000

1.8 mm Blade Lancets (Box of 100)

These spring loaded lancets are ideal for drawing blood for use with the STAR Tubes. The green 1.8 mm version is designed for pediatric use.

420001

2.3 mm Blade Lancets (Box of 100)

The blue 2.3 mm Blade Lancets are designed for use on adult patients.

429580

Printer Paper (3/pk)

Three rolls of 2.5 inch thermal paper for use with the STAR's internal printer.

429581

Fuses, 250 V 3.15 A 5x20 mm (5/pk)

Package of five replacement fuses for the STAR.

NOTE: For use with STAR analyzers manufactured on or before March 31, 2013

429640

Patient ID Upgrade Keyboard

Used to enter Patient and Operator ID numbers into the STAR.

NOTE: For use with STAR analyzers manufactured on or before March 31, 2013

429645

Replacement Patient ID Keyboard

Replacement keyboard used to enter Patient and Operator ID numbers into the STAR.

NOTE: For use with STAR analyzers manufactured after April 1, 2013

429642

Patient ID Barcode Scanner

Used for scanning barcodes found on the box of QBC STAR Tubes.

NOTE: For use with STAR analyzers manufactured after April 1, 2013

421634**QBC Line Cord (USA)**

Replacement line cord for STAR analyzers in use within the United States.

421554**QBC Line Cord (U.K.)**

Replacement line cord for STAR analyzers in use within the United Kingdom.

421551**QBC Line Cord (Europe)**

Replacement line cord for STAR analyzers in use within Europe.

9000-900-000**Carrying Case (Small)**

This hard plastic case comes with wheels and foam inserts to easily transport and protect the STAR during travel. Dimensions - 20.5 inches long x 20.5 inches high x 19.25 inches deep. Weight - 30 pounds.

429627**Carrying Case (Large)**

Protect the STAR during transportation with this hard shell waterproof case. Includes four wheels and convenient handle for moving the case. Large enough to accommodate the STAR, external power supply, external keyboard, Operators manual and other accessories used in the operation of the STAR. Dimensions – 31.28 inches long x 24.21 inches high x 17.48 inches deep. Weight – 36 pounds

429644**Replacement External Power Supply**

Replacement power supply for use with STAR analyzers manufactured after April 1, 2013.

6000-200-149**USB 8GB Memory Stick**

Used for downloading support file.

