

## THE DRUCKER ADVANTAGE



### A TRUSTED BRAND

More than 75% of blood samples collected in physicians' offices and patient centers are processed in Drucker centrifuges. Drucker Diagnostics is the exclusive supplier to large reference laboratories. Drucker enjoys active partnerships with all major brands of blood collection tubes.

The health of astronauts on the NASA space station is currently monitored with a Drucker Diagnostics centrifuge.



### PRE-ANALYTIC PRODUCTS FOR THE CLINICAL LAB

Drucker offers a broad line of pre-analytic products for the clinical lab that are backed by a superb team of professionals with clinical lab expertise.

Drucker products support the best pre-analytic practices for sample quality with horizontal centrifugation and protection from extreme temperatures during storage and shipment as recommended by CLSI.



### EXCELLENT SERVICE AND SUPPORT

Drucker offers a free 30-day evaluation on its products (U.S. only) without purchase commitment. Drucker ships directly from its U.S. Manufacturing facilities. Drucker Diagnostics is known for the excellence of its customer service and technical support which are both based in the US.



### AFFORDABLE, QUALITY PRODUCTS

Drucker products are made in the USA and backed by a 2-year warranty for the entire product that is extendable to 5 years, plus a standard lifetime rotor warranty.

Drucker is an ISO 13485 certified company.

## A CENTRIFUGE LINE TO SIMPLIFY AND STREAMLINE YOUR LAB OPERATION

Achieve operating efficiency and standardization by leveraging the consistency and value of Drucker Diagnostics products across the IDN. Empower the lab to replace the existing assortment of legacy centrifuges of different brands and different models with Drucker centrifuges. Drucker Diagnostics offers versatile units for routine, STAT, outreach, and mobile use to provide a solution for every location.

Designed with simplicity and ease of use in mind, Drucker Diagnostics pre-analytic products streamline training and adoption with consistent interfaces across the product line. Further increase operational efficiency with a single point of contact for tech support or manufacturer service and a streamlined set of spare parts. All these factors combine to deliver a significant cost and time savings.

Only one company offers a comprehensive and cohesive pre-analytical product collection to serve diagnostic needs across the IDN. Drucker would be happy to create a customized package for your network's unique needs.

✓ EASY TO CHOOSE ✓ EASY TO USE ✓ BUILT TO LAST ✓ A GREAT VALUE



druckerdiagnostics.com  
sales@druckerdiagnostics.com

+1-866-265-1486 (U.S. only)  
+1-814-692-7661



Form No. 1005 Rev. A

# Drucker Diagnostics

# Centrifuges & Lockboxes





## DASH for STAT

## HORIZON for ROUTINE

## OUTREACH and MOBILE

## FIXED ANGLE

## SELECT BY TUBE SIZE

### ALL BLOOD TUBES 75-100 mm, up to 10 mL

Horizontal centrifugation, reinforced carbon fiber components  
Lid lighting indicates cycle status (ready, running, done)

#### SET & LOCK

3 Settings (3 min, 5 min, Custom)  
Simple 2-Button Operation

6 Tubes



DASH Apex 6 | 4000 xg Max

12 Tubes



DASH Apex 12 | 4000 xg Max

24 Tubes



DASH Apex 24 | 3000 xg Max

#### FLEX

Digital Display  
10 Programmable Settings



DASH Flex 12 | 4000 xg Max

### ALL BLOOD & URINE TUBES 75-100 mm, up to 10 mL

Horizontal centrifugation, engineered composite components  
Lid lighting indicates cycle status (ready, running, done)

#### SET & LOCK

3 Settings  
Simple 2-Button Operation

6 Tubes



HORIZON 6 | 2000 xg Max

12 Tubes



HORIZON 12 | 2000 xg Max

24 Tubes



HORIZON 24 | 2000 xg Max

#### FLEX

Digital Display  
10 Programmable Settings

6 Tubes



HORIZON 6 Flex | 2000 xg Max

6 Tubes



HORIZON 12 Flex | 2000 xg Max



HORIZON 24 Flex | 2000 xg Max

### ALL BLOOD TUBES 75-100 mm, up to 10 mL

Horizontal centrifugation, single speed, designed to preserve separated specimen integrity during travel

#### OUTREACH

Outreach-Physician Offices  
(Chemistry & Coag)

6 Tubes



Model 642E | 1600 xg Max

6 Tubes



Model 642B | 1600 xg Max



SampleSafe™ Front-Load  
Insulated Lockbox

#### MOBILE

Mobile Phlebotomy  
(Chemistry)



Model 642M | 1600 xg Max



SampleSafe™ Top-Load  
Insulated Lockbox

### BLOOD & URINE TUBES 75-100 mm, up to 125 mm

Fixed angle rotor for coag, PRP, and urine; composite components  
Lid lighting indicates cycle status (ready, running, done)

#### SET & LOCK

3 Settings  
Simple 2-Button Operation

6 Tubes



HORIZON 6 FA | 1850 xg Max

#### FLEX

Digital Display  
10 Programmable Settings



HORIZON 6 Flex FA | 1850 xg Max

#### ANALOG

6 Tubes



Model 614B | 1100 xg Max  
For Blood Only

### BLOOD TUBES

TUBE SIZE	HORIZONTAL	FIXED ANGLE
75-100 mm up to 10 mL	All Horizontal Centrifuges	HORIZON 6 Flex FA, HORIZON 6 FA, 614B
125 mm up to 15 mL	HORIZON 24, HORIZON 24 Flex	HORIZON 6 Flex FA, HORIZON 6 FA, 614B
6 x 50 mL (with alternate carrier)	HORIZON 24, HORIZON 24 Flex	HORIZON 12* or 24*, HORIZON 12 Flex* or 24 Flex*
Pediatric (with extenders)	DASH Apex, DASH Flex 12	N/A

\* with alternate rotor

### URINE TUBES

TUBE SIZE	HORIZONTAL	FIXED ANGLE
75-100 mm up to 10 mL	All HORIZON Centrifuges	HORIZON 6 Flex FA, HORIZON 6 FA
125 mm up to 15 mL	HORIZON 24, HORIZON 24 Flex	HORIZON 6 Flex FA, HORIZON 6 FA

All centrifuges available to evaluate at no cost for 30 days (U.S. only)

Centrifuges arrive ready-to-run with rotor and tube holders included