

fisherbrand

DASH

DASH FLEX 12 & FISHERBRAND 12 STAT CENTRIFUGES

Always know when your STAT samples are ready: the unique flashing lid lights (patent pending) make it easy to track cycle status and keep TAT down. With your choice of digital display or repeatable two-button operation, both units can be programmed with your validated cycles and arrive preset with common STAT spin cycles for convenience.



Easy to Monitor

LED lid lighting indicates when the centrifuge is ready (off), running (on), or done (flashing), reducing forgotten samples and keeping your average TAT down. The DASH Flex 12's digital display shows cycle parameters and time remaining.

Easy to Use

Program up to 10 custom time, RPM, and g-force settings with the DASH Flex 12. Or, set and lock your three preferred cycles using the Fisherbrand 12's two-button interface. Both arrive pre-programmed with the most common STAT settings and can spin 75-100 mm (3-10 mL) tubes with the same rotor and universal tube holders.

Eliminate Wait Time

DASH centrifuges are so compact that they can be used in an array in accessioning. The array will allow spinning samples as soon as they arrive, further shortening your TAT. For a more compact unit, consider the 6 tube capacity DASH Apex 6.

Safe and Durable

Enjoy years of worry-free operation with the DASH's brushless motor and rugged, reinforced carbon fiber components. Enhance lab safety with imbalance protection and a clear, shatterproof lid, which locks until the centrifuge has stopped completely.

A Great Value

Each DASH centrifuge comes fully equipped and ready to run with a horizontal rotor, universal tube holders, and 2-year warranty included in the affordable price.



Choose the DASH Flex 12 for flexibility & cycle display. Choose the Fisherbrand 12 for repeatability with set and lock controls.

Tubes	Applications	Centrifugation	Included	Warranty
12 75-100 mm 3-10 mL	STAT Chemistry, Coag, PPP	Horizontal 4,000 xg Max RCF 5,200 Max RPM	Rotor, Buckets or Tube Holders	2 years Rotor: lifetime warranty Extended warranty available

Designed, built, and supported in the USA

DASH FLEX 12 & FISHERBRAND 12



Lid lights indicate cycle status: ready, running, done

LED lights clearly indicate selected cycle

Program and name your custom validated cycles
Store up to 10 settings



View RPM and remaining spin time on the digital display

Select a preset cycle, program your validated cycles, or lock to a single setting

Free 30-Day Evaluation (U.S. Only)

Specifications	DASH Flex 12	Fisherbrand 12
Sound level	61 dBA	61 dBA
Dimensions (W x D x H)	13 x 15 x 9 in (33 x 38 x 23 cm)	13 x 15 x 9 in (33 x 38 x 23 cm)
Weight	30 lbs (13.7 kg)	30 lbs (13.7 kg)
Power	95 to 253 VAC – 50/60 Hz	
Part Numbers	DASH Flex 12	Fisherbrand 12
Centrifuge	FCAT # 22029678	FCAT # 22029557A
Replacement 75 & 100 mm black tube holder	03-1-0007-0046	03-1-0007-0046
Tube Adapters, 0.5 to 1 mL	FCAT # 22029426 (1 qty)	FCAT # 22029426 (1 qty)
Tube Adapters, 1.5 to 2 mL	FCAT # 22029427 (1 qty)	FCAT # 22029427 (1 qty)
1 yr Warranty Extension	03-0-0007-0003	03-0-0007-0003
2 yr Warranty Extension	03-0-0007-0004	03-0-0007-0004
3 yr Warranty Extension	03-0-0007-0005	03-0-0007-0005
Custom Settings *	00-100-100-101	00-100-100-101

* Add this part number to your centrifuge order to have your custom validated cycles pre-programmed at the factory

Factory Settings	DASH Flex 12			Fisherbrand 12		
	RPM	Time	G-Force	RPM	Time	G-Force
1st Setting/Cycle	5,200	3 min	4,000 xg	5,200	3 min	4,000 xg
2nd Setting/Cycle	5,200	5 min	4,000 xg	5,200	5 min	4,000 xg
3rd Setting/Cycle	3,600	7 min	2,000 xg	3,600	7 min	2,000 xg

FDA LISTED



RoHS
Compliant



ISO13485
certified



druckerdiagnostics.com
sales@druckerdiagnostics.com

+1-866-265-1486 (U.S. only)
+1-814-692-7661

