













Drucker Centrifuge Comparison Chart for Fisher Horizontal Routine



Model	6-Tube Routine Centrifuges				12-Tube Routine Centrifuges		24-Tube Routine Centrifuges	
	642E	642B*	HORIZON 6	HORIZON 6 Flex**	HORIZON 12	HORIZON 12 Flex**	HORIZON 24	HORIZON 24 Flex**
Application	Routine	Routine	Routine	Routine	Routine	Routine	Routine	Routine
Rotor Type	Horizontal	Horizontal	Horizontal	Horizontal	Horizontal	Horizontal	Horizontal	Horizontal
Urine			5 mins	5 mins	5 mins	5 mins	5 mins	5 mins
Sodium Citrate (Coag) PPP	15 mins	15 mins	15 mins	15 mins	15 mins	15 mins	15 mins	15 mins
Plasma Serum (Chemistry)	10 mins	10 mins	10 mins	10 mins	10 mins	10 mins	10 mins	10 mins
Digital Display				✓		✓		✓
Preset Time	10 mins	n/a	10/15/5 mins	n/a	10/15/5 mins	n/a	10/15/5 mins	n/a
Preset or Max RPM	set at 3400	set at 3400	3600/3300/1900	3800 max	3500/3200/1800	3700 max	3300/3000/1700	3400 max
Preset or Max G-Force RCF	set at 1600	set at 1600	1800/1500/500	2000 max	1800/1500/500	2000 max	1800/1500/500	2000 max
Adjustable Speed				✓		✓		✓
Adjustable Time		✓		✓		✓		✓
Up to 17 x 75 or 100 mm or 10 mL	6 tube	6 tube	6 tube	6 tube	12 tube	12 tube	24 tube	24 tube
Up to 17 x 125 mm or 15 mL							12 tube	12 tube
50 mL							6 tube	6 tube
Voltage	115 VAC	115 VAC	95-253 VAC	95-253 VAC	95-253 VAC	95-253 VAC	95-253 VAC	95-253 VAC
Footprint	1.00 sq. ft	1.00 sq. ft	1.00 sq. ft	1.00 sq. ft	1.35 sq. ft	1.35 sq. ft	1.7 sq. ft.	1.7 sq. ft.
Cat. No.	22029375	22029379	22029670	22029672	22029674	22029676	22029675	22029677



STAT, Mobile, Fixed Angle Routine



Model	6, 12, and 24 Tube STAT Centrifuges			
	DASH Apex 6	Fisherbrand 12	DASH Flex 12**	DASH Apex 24
Application	STAT	STAT	STAT	STAT
Rotor Type	Horizontal	Horizontal	Horizontal	Horizontal
Urine				
Sodium Citrate (Coag) PPP	5 mins	5 mins	5 mins	6 mins
Plasma Serum (Chemistry)	7 mins	7 mins	7 mins	7 mins
Digital Display			✓	
Preset Time	3/5/7 mins	3/5/7 mins	n/a	5/6/7 mins
Preset or Max RPM	5300/5300/3800	5200/5200/3600	5200 max	4200/4200/3500
Preset or Max G-Force RCF	4000/4000/2000	4000/4000/2000	4000 max	3000/3000/2000
Adjustable Speed			✓	
Adjustable Time			✓	
Up to 17 x 75 or 100 mm or 10 mL	6 tube	12 tube	12 tube	24 tube
Up to 17 x 125 mm or 15 mL				
50 mL				
Voltage	95-253 VAC	95-253 VAC	95-253 VAC	95-253 VAC
Footprint	1.00 sq. ft	1.35 sq. ft	1.35 sq. ft	1.7 sq. ft.
Cat. No.	22029440	22029557A	22029678	22029445

Mobile Centrifuge	
Model	642M
Application	Mobile
Rotor Type	Horizontal
Urine	
Sodium Citrate (Coag) PPP	15 mins
Plasma Serum (Chemistry)	10 mins
Digital Display	
Preset Time	10 mins
Preset or Max RPM	set at 3400
Preset or Max G-Force RCF	set at 1600
Adjustable Speed	
Adjustable Time	
Up to 17 x 75 or 100 mm or 10 mL	6 tube
Up to 17 x 125 mm or 15 mL	
50 mL	
Voltage	12V DC
Footprint	1.00 sq. ft
Cat. No.	22029587

Fixed Angle Routine Centrifuges			
Model	HORIZON 6 FA	HORIZON 6 Flex FA**	614B
Application	Routine	Routine	Routine
Rotor Type	Fixed Angle	Fixed Angle	Fixed Angle
Urine	5 mins	5 mins	
Sodium Citrate (Coag) PPP	15 mins	15 mins	
Plasma Serum (Chemistry)	10 mins	10 mins	15 mins
Digital Display		✓	
Preset Time	10/15/5 mins	n/a	n/a
Preset or Max RPM	3900/3800/2200	3900 max	set at 3250
Preset or Max G-Force RCF	1600/1500/500	1850 max	set at 1100
Adjustable Speed		✓	
Adjustable Time		✓	✓
Up to 17 x 75 or 100 mm or 10 mL	6 tube	6 tube	6 tube
Up to 17 x 125 mm or 15 mL	6 tube	6 tube	6 tube
50 mL			
Voltage	95-253 VAC	95-253 VAC	115 VAC
Footprint	1.00 sq. ft	1.00 sq. ft	1.00 sq. ft
Cat. No.	22029671	22029673	22029380

