

# Front-Load Lockbox



## Can pay for itself in one day of bad weather

Heat and cold can both render specimens unusable. If a lab processes at least six samples per day, one day of extreme temperatures without the SampleSafe lockbox could ruin all six samples, at a cost of up to \$350 in re-draws. With a SampleSafe protecting the specimens, no re-draws are required - a \$350 savings, more than double the retail price of the lockbox.

### Time to Overheat Comparison

100°F (38°C) Outside Temperature, with Cold Pack



Insulated Steel Box



SampleSafe™ Lockbox

### Time to Freezing Comparison

15°F (-9°C) Outside Temperature



Insulated Steel Box



SampleSafe™ Lockbox

## FRONT LOAD LOCKBOX

Height ..... 14 in. / 35.56 cm  
 Width ..... 15 in. / 38.1 cm  
 Depth ..... 8 in. / 20.32 cm  
 Weight ..... 11.5 lbs. / 5.22 kg  
 Internal Volume ..... 570 in<sup>3</sup> / 9,340.62 cm<sup>3</sup>  
 Warranty ..... 5 years



**Front Load Lockbox**


P/N: 00-078-009-000

Also available:



**Top Load Lockbox**








P/N: 00-071-009-000



Secure anti-theft mounting options (door hanging, under door, or security cable).

Lid can be customized with your logo.

## OPTIONAL ACCESSORIES

ANTI-THEFT DEVICES		PART NUMBER
	Door Hanging Bracket - 24"	02-004-0-0007
	Door Hanging Bracket - 36"	00-078-109-000
	Under Door Bracket Kit	02-004-0-0002
	Anchor Bracket Kit (Use with Cable)	02-004-0-0003
	Anchor Cable (6 Foot)	03-1-0012-0054
	Wall Mounting Kit (Bracket & Screws)	02-004-0-0009
FOR ADDITIONAL PERFORMANCE		PART NUMBER
	Cold Pack (6" x 9" x .8")	03-1-0012-0053