

# SERO 12

## SERO 12 VS SERO-FUGE *BLOOD BANKING CENTRIFUGES*


**SERO 12**

**Sero-Fuge™**

Lid Light Flashes on Cycle Completion*	✓	
Accepts Primary Tube Sizes without Inserts	✓	✓
Recommended by Tube Manufacturers	✓	
Available Extended Warranty	✓	
Warranty	2 yrs	1 yr
Weight (lbs)	27	28
Footprint (sq ft)	1.35	1.53
Number of Tubes	6, 12	6, 12
Numbered Tube Positions	✓	✓
Digital Display	✓	✓
User Programmable Fixed Cycles	10	0
Maximum Speed	3,600 RPM	3,600 RPM
Maximum RCF	1,040 xg	1,040 xg

\* Patent Pending

## Uses All **BD SeroFuge** Rotors

