

BOOST 2+ FLEX

RECOMMENDED BY VETSTEM

Spin large tubes in a small centrifuge with the BOOST 2+ Flex. Fitted with our largest buckets to hold tubes up to 1.9"W x 4.5"H (4.8 x 13.7 cm). A lightweight fit for any space with a small footprint of 1.35 ft² (1,254 cm²).



Custom Inserts

Two inserts (pictured right) are included with each purchase.



Fully Customizable

The BOOST 2+ Flex offers a variety of extensive customization options. Personalize the front panel label and LED lid lights to match your unique branding, or pre-program custom RPM and g-force settings to meet your preferred processing specifications.

Spins Your Cycles

Program and recall up to 10 preset cycles with the BOOST 2+ Flex's backlit digital display. Can be set to run preset cycles only, or may be locked to just one preset cycle if desired.

Fits Your Tubes

Experience unmatched convenience with BOOST 2+ Flex's unique insert system. Simply select the correct inserts for your tube sizes and switch them out as needed to spin every tube.

Safe and Durable

Enjoy years of worry-free operation with BOOST's imbalance protection and shatterproof lid, which can only be opened when the rotor is stopped - even if the unit loses power. Rigorous electrical testing ensures our centrifuges meet hospital equipment safety standards.

A Great Value

Unbeatable quality and flexibility for any budget. Each BOOST 2+ Flex centrifuge is ready to run with swing-out horizontal rotor, buckets, inserts, and 2-year warranty included in the affordable price.

Applications	Max Speed	Max Capacity	Included	Warranty
Variable Tube Sizes Custom Disposables	4,400 RPM 3,000 xg RCF	2 buckets Maximum load 140.0 grams/ bucket including inserts	Horizontal Rotor, Tube Holders, Inserts	2 years Rotor: lifetime warranty Extended warranty available

BOOST 2+ Flex



- LED lid lights indicate cycle status: ready, running, done
- Easily spin a wide variety of tubes with swappable inserts
- View RPM, RCF, and remaining spin time on the digital display
- Program and name your custom validated cycles
- Store up to 10 settings
- Lightweight and compact for an easy fit in any clinic
- QR code for easy access to support videos, manuals, and more

Part Numbers	
Centrifuge & 2 Inserts	00-783-248-000
Extra Inserts (2 qty)	03-1-0007-0118K

Specifications	
Sound Level	70 dBA
Max Speed/G-Force	4,400 RPM/3,000 xg RCF
Dimensions (W x D x H)	13 x 15 x 9 in (33 x 38 x 23 cm)
Weight	34 lbs (15.4 kg)
Power	95 to 253 V – 50/60 Hz
Resting Tube Angle	20° (vertical)
Tube Spin Angle	90° (horizontal)
Maximum Load	140g per bucket



Scan to Purchase

Drucker Diagnostics centrifuges use medical grade components and are tested to rigorous international standards. These certifications allow our products to be used in hospitals and medical facilities around the world.



druckerdiagnostics.com
sales@druckerdiagnostics.com

+1-866-265-1486 (U.S. only)
+1-814-692-7661

