




# Experience the world's best-selling clinical laboratory centrifuge



## Chemistry, coag, and PPP

The 642E is the ideal choice for doctor's offices, patient service centers, and outreach locations.




## 1,600 xg 3,400 RPM

The 642E features a preset runtime of 10 minutes at 1,600 xg.



## 6-tube capacity

Spin 75 & 100 mm (3-10 mL) tubes in the included tube holders.



## Horizontal centrifugation

The 642E's patented swing-out horizontal rotor design produces high-quality samples.

Please contact your [Brand] representative for more information at ###-###-####.