



# Produce platelet poor plasma in just 3 minutes




00-676-260-000

✓ Set-and-lock controls


✓ Durable polymer components

✓ Affordable price




## Streamline coag processing

The DASH Coag is ideal for hospitals that handle high volumes of platelet poor plasma (PPP)/coag testing.




## 4,400 xg 6,600 RPM

The DASH Coag features three carefully selected presets for coagulation.



## 12-tube capacity

Spin tubes up to 75 mm (4 mL) with the included rotor and tube holders.



The **fixed angle rotor** remains at a 45° angle for optimal PPP formation.

Please contact your Boekel representative for more information



PROUDLY MADE IN THE USA