



CA676104


# Produce platelet poor plasma in just 3 minutes



✓ Set-and-lock controls


✓ Durable polymer components

✓ Affordable price




## Streamline coag processing

The Cardinal Health Brand™ DASH Coag is ideal for hospitals that handle high volumes of platelet poor plasma (PPP)/coag testing.




4,400 xg  
6,600 RPM

The Cardinal Health Brand DASH Coag features three carefully selected presets for coagulation.



## 12-tube capacity

Spin tubes up to 75 mm (4 mL) with the included rotor and tube holders.



The fixed angle rotor remains at a 45° angle for optimal PPP formation.

[cardinalhealth.com](http://cardinalhealth.com)

© 2023 Cardinal Health. All Rights Reserved. CARDINAL HEALTH and the Cardinal Health LOGO are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. (09/2023)

