



76795-588


# Produce platelet poor plasma in just 3 minutes



✓ Set-and-lock controls


✓ Durable polymer components

✓ Affordable price




## Streamline coag processing

The VWR® STAT Coag Clinical Centrifuge is ideal for hospitals that handle high volumes of platelet poor plasma (PPP)/coag testing.




## 4,400 xg 6,600 RPM

The VWR® STAT Coag Clinical Centrifuge features three carefully selected presets for coagulation.



## 12-tube capacity

Spin tubes up to 75 mm (4 mL) with the included rotor and tube holders.



The **fixed angle rotor** remains at a 45° angle for optimal PPP formation.

Please contact your VWR representative for more information.