



DASH Centrifuges - High Performance STAT Processing

| Model | Rotor Type | Interface | Max Speed / RCF | Up to 16 x 75/ 100 mm or 10 mL | Up to 16 x 125 mm or 15 mL | 50 mL | Part Number |
|--------------|-------------|---------------------|-----------------|--------------------------------------|----------------------------------|-------|----------------|
| DASH Apex 4 | Fixed Angle | 3-Setting | 8500 / 4440 | 4 (75 mm only) | - | - | 00-053-009-000 |
| DASH Coag | Fixed Angle | 3-Setting | 6600 / 4400 | 12 | - | - | 00-676-009-001 |
| DASH Apex 6 | Horizontal | 3-Setting | 5300 / 4000 | 6 | - | - | 00-076-009-004 |
| DASH Flex 6 | Horizontal | 10-Setting, Digital | 5300 / 4000 | 6 | - | - | 00-176-009-000 |
| DASH Apex 12 | Horizontal | 3-Setting | 5200 / 4000 | 12 | - | - | 00-083-009-003 |
| DASH Flex 12 | Horizontal | 10-Setting, Digital | 5200 / 4000 | 12 | - | - | 00-183-009-000 |
| DASH Apex 24 | Horizontal | 3-Setting | 4300 / 3150 | 24 | 12 | 6 | 00-084-009-005 |
| DASH Flex 24 | Horizontal | 10-Setting, Digital | 4500 / 3450 | 24 | 12 | 6 | 00-184-009-000 |

SERO 12 - Blood Banking Centrifuge A drop-in replacement for the BD SeroFuge

| Model | Rotor Type | Interface | Max Speed / RCF | 12 Place Rotor Up to 75mm | 6 Place Rotor Up to 100 mm | Part Number |
|---------|-------------|---------------------|-----------------|------------------------------|-------------------------------|----------------|
| SERO 12 | Fixed Angle | 10-Setting, Digital | 3600 / 1040 | 6 capped 12 uncapped | 6 | 00-483-009-001 |

HORIZON Centrifuges Quality Routine Chemistry

| Model | Rotor Type | Interface | Max Speed / RCF | Up to 16 x 75/ 100 mm or 10 mL | Up to 16 x 125 mm or 15 mL | 50 mL | Part Number |
|-------------------|-------------|---------------------|-----------------|--------------------------------------|----------------------------------|-------|----------------|
| HORIZON 6 | Horizontal | 3-Setting | 3800 / 2000 | 6 | - | - | 00-276-009-000 |
| HORIZON 6 FA | Fixed Angle | 3-Setting | 3900 / 1850 | 6 | 6 | - | 00-276-009-001 |
| HORIZON 6 Flex | Horizontal | 10-Setting, Digital | 3800 / 2000 | 6 | - | - | 00-376-009-000 |
| HORIZON 6 Flex FA | Fixed Angle | 10-Setting, Digital | 3900 / 1850 | 6 | 6 | - | 00-376-009-001 |
| HORIZON 12 | Horizontal | 3-Setting | 3700 / 2000 | 12 | - | - | 00-238-009-000 |
| HORIZON 12 Flex | Horizontal | 10-Setting, Digital | 3700 / 2000 | 12 | - | - | 00-383-009-000 |
| HORIZON 24 | Horizontal | 3-Setting | 3400 / 2000 | 24 | 12 | 6 | 00-284-009-000 |
| HORIZON 24 Flex | Horizontal | 10-Setting, Digital | 3400 / 2000 | 24 | 12 | 6 | 00-384-009-000 |

600 Series Centrifuges Compact Routine Chemistry

| Model | Rotor Type | Interface | Max Speed / RCF | Up to 16 x 75/ 100 mm or 10 mL | Up to 16 x 125 mm or 15 mL | 50 mL | Part Number |
|-------|-------------|-----------------|-----------------|--------------------------------------|----------------------------------|-------|----------------|
| 614B | Fixed Angle | 1-Setting, Dial | 3250 / 1200 | 6 | 6 | - | 00-004-009-001 |
| 642B | Horizontal | 1-Setting, Dial | 3380 / 1600 | 6 | - | - | 00-016-009-001 |
| 642E | Horizontal | 1-Setting | 3380 / 1600 | 6 | - | - | 00-018-009-001 |
| 642M | Horizontal | 1-Setting | 3400 / 1650 | 6 | - | - | 00-077-009-000 |