



TrueBond Duet

Compact Dual Rotor Centrifuge

Part # 00-056-009-000

The TrueBond Duet is Drucker Diagnostics' most compact animal health centrifuge, designed to spin microhematocrit tubes, blood, and urine—all in one unit. Its removable rotors allow you to handle a wide variety of veterinary samples with ease.

Easy to Use

The simple two-button interface makes it easy to select a cycle. Use all five presets for routine veterinary testing, or lock to a single cycle for maximum repeatability.

Compact Size

Weighing only 4.3 lbs (1.95 kg) and measuring 7 in x 7 in x 7 in (18 cm x 18 cm x 18 cm), the TrueBond Duet is our smallest veterinary centrifuge.

Safe and Durable

The TrueBond Duet's auto-pop lid locks securely during a run, and pops open automatically when the cycle is done. The brushless motor, reinforced carbon fiber components, and two-year calibration guarantee ensure years of operation with no routine maintenance.

A Great Value

The TrueBond Duet includes an easily swappable fixed angle rotor, microhematocrit rotor, microhematocrit reader card, and a 2-year warranty in the affordable price.



Removable Rotors

The TrueBond Duet includes two quick swap rotors, allowing a seamless transition between a 2-place fixed angle rotor (left) for 1.5-2 mL blood or urine tubes, to a 12-place 40 mm microhematocrit rotor (right).

Applications	Max Speed	Max Capacity	Included	Warranty
Microhematocrit, Urine, Blood	15,800 Max RPM 12,000 xg Max RCF (1.5 - 2 mL) 13,700 xg Max RCF (Microhematocrit)	2 x 1.5 - 2 mL Blood or Urine Tubes 12 x 40mm Microhematocrit Tubes	2-Place Fixed Angle Rotor, 12-Place Microhematocrit Rotor Microhematocrit Reader Card	2 years Rotor: lifetime warranty Extended warranty options available

Designed, built, and supported in the USA



5 default presets for routine veterinary testing

QR code for easy access to support manuals and more

Select presets using the 2-button interface

Press to select cycle

Free 30-Day Evaluation (U.S. Only)



Standard Preset Cycles		Fixed Angle Rotor (1.5 - 2 mL)	MHCT Rotor (Microhematocrit)
Quick	30 sec @ 15800 RPM	12000 xg	N/A
Normal	1.5 min @ 15800 RPM	12000 xg	N/A
Long	2.5 min @ 15800 RPM	12000 xg	N/A
Urine	45 sec @ 9800 RPM	3900 xg	N/A
Microhematocrit	2 min @ 15800 RPM	N/A	13700 xg

Part Numbers		
Included	Centrifuge	00-056-009-000
	2-Place Fixed Angle Rotor	02-001-0-0029
	Microhematocrit Rotor	02-001-0-0030
	Veterinary Microhematocrit Reader Card	03-1-0012-0232
Optional	1 Year Extended Warranty	03-0-0007-0001
	2 Year Extended Warranty	03-0-0007-0002
	3 Year Extended Warranty	03-0-0007-0006

Specifications	
Max RPM	15800 RPM
Max G-Force (1.5 - 2 mL)	12000 xg
Max G-Force (Microhematocrit)	13700 xg
Dimensions (W x D x H)	7 x 7 x 7 in 18 x 18 x 18 cm
Weight	4.3 lbs (1.95 kg)
Power	100 - 240 (+/- 10%) VAC - 50/60 Hz



druckerdiagnostics.com
sales@druckerdiagnostics.com

+1-866-265-1486 (U.S. only)
+1-814-692-7661

Protected by U.S. Patents #6.811.531, #D718,463S #D734,489S Other Patents Pending.