

**Drucker Diagnostics:** a trusted centrifugation leader manufacturing products with the quality and value your lab deserves. All Drucker Diagnostics centrifuges are designed, built, and supported in the USA and include a 2-year warranty.

**DASH STAT Series**

Built to spin STAT, DASH centrifuges deliver quick spins when it matters most.



DASH Apex 4

DASH Coag

DASH Flex 6\*

DASH Flex 12\*

DASH Flex 24\*

*\*Also available with two-button interface*

**HORIZON Routine Series**

Quality and value for routine chemistry.



HORIZON 6 Flex\*

HORIZON 6 Flex FA\*

HORIZON 12 Flex\*

HORIZON 24 Flex\*

*\*Also available with two-button interface*

**Collection Site Series**

Compact and durable, designed for years of routine chemistry use.



Model 642E

Model 642M

Model 642B

Model 614B

**Blood Banking**


A drop-in replacement for the BD Sero-Fuge™ that uses the exact same rotors.



SERO 12



All centrifuges are available for a free 30-day evaluation. Visit: [druckerdiagnostics.com/evaluation-program](http://druckerdiagnostics.com/evaluation-program)



## DASH Centrifuges - High Performance STAT Processing

Model	Rotor Type	Interface	Max Speed RPM / RCF	Up to 13 x 75 mm or 5 mL	Up to 16 x 100 mm or 10 mL	Up to 16 x 125 mm or 15 mL	50 mL	Part Number
DASH Apex 4	Fixed Angle	2-Button, 3-Setting	8500 / 4440	4	-	-	-	00-053-009-000
DASH Coag	Fixed Angle	2-Button, 3-Setting	6600 / 4400	-	12	-	-	00-676-009-001
DASH Apex 6	Horizontal	2-Button, 3-Setting	5300 / 4000	-	6	-	-	00-076-009-004
DASH Flex 6	Horizontal	Digital, 10-Setting	5400 / 4150	-	6	-	-	00-176-009-000
DASH Apex 12	Horizontal	2-Button, 3-Setting	5200 / 4000	-	12	-	-	00-083-009-003
DASH Flex 12	Horizontal	Digital, 10-Setting	5200 / 4000	-	12	-	-	00-183-009-000
DASH Apex 24	Horizontal	2-Button, 3-Setting	4300 / 3150	-	24	12	6	00-084-009-005
DASH Flex 24	Horizontal	Digital, 10-Setting	4500 / 3450	-	24	12	6	00-184-009-000

(adapters available for a variety of tube sizes)

## SERO 12 - Blood Banking Centrifuge A drop-in replacement for the BD SeroFuge

Model	Rotor Type	Interface	Max Speed RPM / RCF	12 Place Rotor Up to 75mm	6 Place Rotor Up to 100 mm	Part Number
SERO 12	Fixed Angle	Digital, 10-Setting	3600 / 1040	6 capped 12 uncapped	6	00-483-009-001

## HORIZON Centrifuges Quality Routine Chemistry

Model	Rotor Type	Interface	Max Speed RPM / RCF	Up to 16 x 100 mm or 10 mL	Up to 16 x 125 mm or 15 mL	50 mL	Part Number
HORIZON 6	Horizontal	2-Button, 3-Setting	3800 / 2000	6	-	-	00-276-009-000
HORIZON 6 FA	Fixed Angle	2-Button, 3-Setting	3900 / 1850	6	6	-	00-276-009-001
HORIZON 6 Flex	Horizontal	Digital, 10-Setting	3800 / 2000	6	-	-	00-376-009-000
HORIZON 6 Flex FA	Fixed Angle	Digital, 10-Setting	3900 / 1850	6	6	-	00-376-009-001
HORIZON 12	Horizontal	2-Button, 3-Setting	3700 / 2000	12	-	-	00-238-009-000
HORIZON 12 Flex	Horizontal	Digital, 10-Setting	3700 / 2000	12	-	-	00-383-009-000
HORIZON 24	Horizontal	2-Button, 3-Setting	3400 / 2000	24	12	6	00-284-009-000
HORIZON 24 Flex	Horizontal	Digital, 10-Setting	3400 / 2000	24	12	6	00-384-009-000

(adapters available for a variety of tube sizes)

## 600 Series Centrifuges Compact Routine Chemistry

Model	Rotor Type	Interface	Max Speed RPM / RCF	Up to 16 x 100 mm or 10 mL	Up to 16 x 125 mm or 15 mL	50 mL	Part Number
614B	Fixed Angle	Dial, 1-Setting	3250 / 1200	6	6	-	00-004-009-001
642B	Horizontal	Dial, 1-Setting	3380 / 1600	6	-	-	00-016-009-001
642E	Horizontal	2-Button, 1-Setting	3380 / 1600	6	-	-	00-018-009-001
642M	Horizontal	2-Button, 1-Setting	3400 / 1650	6	-	-	00-077-009-000

(adapters available for a variety of tube sizes)



Designed, built, and supported in the USA

+1-866-265-1486 (U.S. only)  
+1-814-692-7661

[druckerdiagnostics.com](http://druckerdiagnostics.com)

email to: [sales@druckerdiagnostics.com](mailto:sales@druckerdiagnostics.com)